

Consult Your Physician Before Using

You should consult with your Doctor/Medical Practitioner if you have any concerns whatsoever before using the Knee Scooter. Neither Knee Scooter Australia or Medical Rehab Pty Ltd nor any employee of these companies is medically qualified to advise you on your suitability to use the Knee Scooter or make an assessment that the Knee Scooter will be appropriate for your purpose or condition.

LIMITED LIFETIME WARRANTY

Knee Scooters Australia Scooters are built to exacting standards and are carefully inspected prior to shipment. Knee Scooters Australia warrants the Knee Scooter frame to be free of defects in materials and workmanship under normal wear and tear from the date of purchase for a Limited Lifetime ("Warranty"). The Warranty extends only to the original purchaser of the Knee Scooter only obligation under the Warranty is to repair or replace the defective product. The warranty does not include labour or shipping charges incurred to repair or replace the defective product.

The Warranty is void if: the product has been repaired or modified by an unauthorised third party; the product is damaged from the use or attachment of unauthorised parts or devices; the product is damaged by the user's failure to adhere to product operating instructions and safety precautions; the product is damaged by negligence, misuse, accident, or improper maintenance; the product is damaged from an act of God or any event outside the control of Knee Scooters Australia

310 Lorimer St Port Melbourne Victoria

Australia 3207

www.kneescooters.com.au

Ph: 0432847576



KC2 / KC2 AT Knee Scooter User Manual



**Please read and understand the instructions carefully before
using the**

" KNEE SCOOTERS AUSTRALIA" SCOOTER

**DO NOT OPERATE THE KNEE SCOOTER WITHOUT FIRST READING
AND UNDERSTANDING THE INFORMATION CONTAINED IN THIS
MANUAL**

GENERAL WARNINGS

Serious bodily injury may occur if the warnings contained in this manual are not strictly adhered to. Failure to follow these safety operating instructions can cause change in stability of the knee scooter increasing the risk of fall or tipping resulting in injury of the user. If you cannot understand the warnings please contact your healthcare provider before use.

DO NOT use this device following the consumption of alcohol or when taking medication that adversely affects coordination and or could impair mental processing

DO NOT attempt to climb or descend slopes greater than 10 degrees.

NEVER attempt to access wet icy or oily surfaces with this device.

DO NOT attempt to operate the knee scooter on public roads, streets or highways. Use only on designated pedestrian walkways and sideways.

AVOID reaching for objects with both hands or extending your reach while utilising this device. Substantial shifts in weight may cause falls or injury.

ALWAYS lock your brake before engaging in hand free use.

DO NOT stand on frame, handlebars or knee rest at anytime.

DO NOT attempt to use this device on rotating escalator. Always use ramps or elevators.

DO NOT use if you have difficulty with balance.

DO NOT use over curbs or uneven ground.

DO NOT use to go faster than a slow walking pace.

DO NOT turn the unit sharply while moving.

DO NOT use in a shower, bath, wet or slippery conditions.

DO NOT use if component of the Knee Scooter is broken or not functioning properly.

PREPARATION FOR USE

It is important that the knee scooter is correctly set up to fit the users height and ability. You should consult your healthcare provider on your limitations and for correct height set up. The knee pad should be set under your knee with your leg bent 90 degrees. Hands should rest on the handle bars with elbows slightly bent.

ASSEMBLY INSTRUCTIONS VIDEO FOUND HERE

<https://youtu.be/eEfSfYr4s1Q>

Riding the knee scooter will provide you with greater mobility than standard crutches. The controls are very similar to a bicycle. with a few important subtle differences. **REMEMBER TO START SLOWLY**, this is a MEDICAL mobility device for the purpose of aiding your recovery.

FORWARD Position your injured leg on the knee pad aligning the kneecap with the front most edge of the pad. Stand erect, with hands comfortably placed on the hand grips. Release the brake and use your free leg to propel your slowly forward. Practice with assistance until you are completely comfortable moving around. Always move at a safe manageable pace.

BACKWARDS Use extreme caution when doing so. Straighten the front wheels and gently push forward with your free leg moving slowly. This is a delicate move that can compromise your balance.

TURNING With both hands resting on the hand grips with correct body position. rotate the handles in the direction desired.

ON A SLOPE (not recommended) Extreme caution must be exercised when traveling down a slope or ramp. Your centre of gravity will change and the possibility of loss of control or injury is greatly increased.

PLEASE LOCK BRAKES WHEN MOUNTING AND DISMOUNTING.