

**NEW  
IMPROVED  
DESIGN**



User Manual

# **Kingston Mobile Chair**

# **COMFORT & MOBILITY**

18005SMG

Updated: February 2023



**Healthcare Solutions**

# Kingston Mobile Chair

The new and improved **Kingston Mobile Chair** is a mobile version from the popular Kingston chair range.

The increased base and backrest frame strength, combined with an adjustable backrest, dual locking castors, and lighter overall weight – now comes with a maximum user weight capacity of 160kg.

The Kingston Mobile Chair provides comfortable, easy-to-clean seating features, with four lockable castors for increased stability and safety.

Armrests are height adjustable and can drop down to enable easy transfer into and out of the chair.

**The Kingston Mobile Chair is not a wheelchair;** it allows users to be moved while seated for short distances in indoor environments.

## Specifications

### 18005SMG

Weight	13kg
Overall Height	915 – 1015mm
Overall Width	600mm
Overall Length	630mm
Seat Width	450mm
Seat Height	460 – 560mm
Seat Depth	440 – 490mm
Max User Weight	160kg
Castors	4 dual castors, swivel locking 75mm / 3"
Materials	Powder coat steel frame, stainless steel arm rests, PU foam seat/backrest textiles, plastic accessories

# Safety Precautions

## Safe Use Environment

Indoor use only. Do not use in wet or slippery areas such as bathrooms or pool areas.


## Maximum Load


The safe working load is 160kg.

## Safety Warnings

1. The product is designed for SHORT TRANSITS (<10m) ONLY
2. Do not use as a wheelchair
3. Not intended to be sat in for extended periods of time
4. Do not use to move items
5. Keep away from wet environments
6. Do not stand on chair
7. Do not sit on armrest
8. Do not apply weight of body onto arm rests
9. Do not adjust legs while a user is seated in the chair



# Components





## Assembly



**Step 01** Insert back rest (A) into frame (B).  
Ensure positioning pins are securely located.



After Assembly

# Assembly


**Step 02** Insert arm rests (A) into frame (B).



After Assembly

## Assembly

**Step 03** Screw arm rest knobs (A) into screw recess holes (B) in armrest outer column. Ensure knobs are tightly fastened.



After Assembly

**Step 04** Once the chair is assembled, adjust legs, arm rests and backrest to desired positions (refer to instructions in this manual).

**Warning:** Do not adjust legs while a user is seated in the chair

# Operating Instructions

## Adjust Chair Height

To adjust the height, press down on the pin on the leg and slide the inner tubing to the preferred height. Repeat the process for each leg until the chair is at the required height.




**Warning:** Check that all safety pins are engaged, failure to do so may result in serious injury.

**Safety pin is not fully engaged**



**Safety pin is fully engaged**



# Operating Instructions

## Adjust Armrest Height

**Step 01** Unscrew the two screw knobs on the armrest you wish to adjust.



**Step 02** To raise the armrest pull it upwards, to lower push it downwards.



**Step 03** When the armrest is at the desired height, tighten the screw knobs.

**Warning:** Always use one hand to hold the armrest that needs to be adjusted to avoid armrest dropping down suddenly when loosening the screw knobs. Do not put your hand under the armrest when adjusting the armrest in case the armrest falls on your hand.


# Operating Instructions

## Lock and Unlock Castors

**Step 01** To lock the castors, push down on the red locking tab on the castor **(A)** until it clicks into place **(B)**.



**Step 02** To unlock the castors, lift the red locking tab by pushing the top lever on the castor until it clicks into place.



# Cleaning and Maintenance

## Inspections

- Ensure adjustable legs are correctly adjusted to match all other legs before every use.
- Ensure adjustable leg positioning pins are securely in place before every use.
- Ensure backrest is correctly adjusted before every use.
- Ensure backrest positioning pins are correctly engaged before every use.
- Ensure arm rests are securely fastened before every use.
- Ensure castors are securely attached before every use.
- Ensure castor brakes are working correctly before every use.

## Cleaning

- Frame may be cleaned with a damp cloth and mild soap.
- Seat, backrest and arm rest pads may be cleaned with a damp cloth and mild soap.
- Do not hose down or use excessive amounts of water to clean chair. Only use a pre-moistened cloth.
- Do not use bleach or harsh chemical agents as these will damage the chair.

# Cleaning and Maintenance

## Servicing

- If you find that your product needs servicing, please contact your dealer for further servicing instructions **OR**
- We recommend a 6-month preventative maintenance service on this product. Please contact our service department or your local dealer to arrange a service programme.

Preventative maintenance service includes:

- Inspection of frame/adjustable legs
- Inspection of seat/backrest pads
- Inspection and clean of castors
- Inspection of arm rest frame/pads
- Inspection of position pins/screw knob handles

## Disposal Instruction

Please dispose of thoughtfully according to local government regulations.

# Warranty Against Defects

Product	Warranty Period
<b>K Care   TUBALCO   R &amp; R</b>	
<b>All Products (Inc. all electrics)</b>	12 Months
<b>Upholstery</b>	12 months
<b>Shower Chairs, Stools &amp; Toilet Seat Raisers (Zinc Treated products)</b>	36 months structural 5 Years anti-corrosion (external tube only)
<b>Shower Chairs, Stools and Toilet Seat Raisers (Aluminium products)</b>	36 months structural 5 Years anti-corrosion (external tube only)
<b>Shower Chairs, Stools and Toilet Seat Raisers (Stainless steel products)</b>	36 months structural 5 years anti-corrosion
<b>KERRY   OXFORD EME   EASYVET</b>	
<b>General Equipment</b>	36 months
<b>KERRY</b>	
<b>Manual Handling Equipment</b>	24 months
<b>Slings</b>	12 months
<b>AIRCOMFORT</b>	
<b>Frame Compact Lift Chair</b>	36 months
<b>Frame Day Beds</b>	36 months
<b>Foam</b>	36 months
<b>Air Bags and Gas Struts</b>	12 months
<b>Castors</b>	12 months
<b>Fabric</b>	12 months
<b>MY COMFORT</b>	
<b>MC100 &amp; MC310 Frames</b>	10 years
<b>MC100 &amp; MC310 Electrics</b>	2 years

1. K Care Healthcare Equipment provides consumers with the following warranty in relation to this Product, in addition to complying with the requirements of any relevant legislation, including the Competition and Consumer Act 2010 (Cth) in Australia and the Consumer Guarantees Act 1993 in New Zealand (the Acts), except where a New Zealand consumer acquires the Product for the purposes of a business.
2. In this warranty, we have used the following definitions:
  - (a) Form means a warranty claim form issued by K Care Healthcare Equipment in respect of Products.
  - (b) K Care Healthcare Equipment or our means K Care Healthcare Equipment,
  - (c) Products means the goods manufactured by K Care Healthcare Equipment (including products manufactured by its contract manufacturers) set out in the table above;
  - (d) Material means a material or component used by K Care Healthcare Equipment in the manufacture of the Products;
  - (e) Retailer means the authorised dealer of Products from whom the Product was purchased;
  - (f) Warranty Period means the warranty period set out in the table above, commencing from the date of purchase of the relevant Product. If the Product or part (as the case may be) is repaired or replaced, there will be no extension to the original warranty period; and
  - (g) Workmanship means the handling, assembly and manufacturing processes performed by or on behalf of K Care Healthcare Equipment in order to manufacture the Products.
3. K Care Healthcare Equipment warrants that for the duration of the Warranty Period, all Products will be free of faults arising from defects in Workmanship or Materials, on the terms and conditions set out in this warranty.
4. K Care Healthcare Equipment undertakes that if during the Warranty Period any Product, or any part of a Product, has failed to operate correctly due to faulty Workmanship or defective Material, it will repair or replace the Product or part (as the case may be) free of charge (excluding freight fees if applicable) provided that the following procedure is met:
  - (a) The consumer must contact the Retailer or K Care Healthcare Equipment upon becoming aware of any defect to a Product. The consumer will then be provided with a Form, which must be completed by the consumer and returned to K Care Healthcare Equipment together with satisfactory proof of purchase.
  - (b) K Care Healthcare Equipment will review the Form to determine whether there is a defect, and if so K Care Healthcare Equipment agrees to (at its option) repair, replace or supply equivalent goods, or pay the cost of any of those remedies to the consumer.
  - (c) If K Care Healthcare Equipment requests the return of the applicable Product or part, the consumer will be responsible for the collection and freight costs of returning that Product or part to K Care Healthcare Equipment.
  - (d) The consumer acknowledges that it is also responsible for the freight costs to deliver any new Product or replacement part to it.
5. To the extent that the supply of Products is a supply to a consumer within the meaning of the applicable Act, the consumer is entitled to consumer guarantees which K Care Healthcare Equipment does not exclude, restrict or modify. In all other respects, and to the extent permitted by law, K Care Healthcare Equipment:
  - (a) Limits its liability for any non-excludable condition or warranty to rectifying any defect at its option, as set out in paragraph 4 (b); and
  - (b) Excludes liability (whether express or implied) of any nature whatsoever for any consequential loss, damage or injury arising as a result of any fault in the Products.
6. This warranty does not extend to damage to Products which occurs during transit or transportation, or which is caused by any abuse, accident or improper installation, connection, use, adjustment or repair or use of goods otherwise than in accordance with instructions issued by K Care Healthcare Equipment.
7. The warranty on Products is waived if any addition or attachment to the Products do not have K Care Healthcare Equipment's approval or are not sold as K Care Healthcare Equipment products. The Products are designed to perform specific tasks under established test loads and unauthorised attachments may produce stresses for which the design is not appropriate.
8. The following applies to consumers who purchased a relevant Product in Australia: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

# Improving your life everyday™

## At Your Service

Sales, administration  
and general enquiries

National 1300 783 783

Your Customer Service Team  
[customerservice@kcare.com.au](mailto:customerservice@kcare.com.au)

**[kcare.com.au](http://kcare.com.au)**

All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information.

