

rehabhire & sales

novis

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

# Procair Prime Advanced

## User Manual



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

# Important notice

Before operating this medical equipment, it is important to read this User Guide and understand the operating instructions and safety precautions. Failure to do so could result in patient injury and/or damage to the product.

We recommend you keep the User Guide near the product.

Therapeutic devices and/or medical equipment should only be used in accordance with manufacturer's instructions and under the consent, supervision and management of a suitably qualified health professional.


If you have any questions, please contact Novis Healthcare on **1300 738 885**.

Novis Healthcare has a policy of continuous product improvement and reserves the right to amend specifications presented in this guide. Information correct at time of production (November 2024).

© 2024 Novis Healthcare. All rights reserved.


## Definitions of symbols Used


The following symbols may appear in this User Guide, on the product, or on its accessories. Some of the symbols represent standards and compliances associated with the control unit and its use.


 Important Information

 Caution


 Electrical Hazard


 Infection Control

 Do Not...

 Double Insulation


 Type BF Applied Part


 Alternating Current


 Manufacturer

 Manufacturing Date

**SN** Serial Number

 Refer to User Guide

 Disposal: Do not dispose of this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.

 IP22 Protection against foreign object and vertically falling water drops.

# Contents

<b>Important Notice</b>	<b>2</b>
<b>Contents</b>	<b>3</b>
<b>System Overview</b>	<b>4</b>
<b>Safety Precautions</b>	<b>6</b>
<b>System Preparation</b>	<b>8</b>
<b>Intended Use</b>	<b>12</b>
<b>Patient Set-Up</b>	<b>13</b>
<b>Operation – Control Handset</b>	<b>14</b>
<b>Operation – Modes of Operation</b>	<b>15</b>
<b>Operation – Quick Twist CPR</b>	<b>15</b>
<b>Operation – Transport Mode</b>	<b>16</b>
<b>Operation – Deflation / Storage</b>	<b>17</b>
<b>Care and Cleaning</b>	<b>18</b>
<b>Troubleshooting</b>	<b>20</b>
<b>Technical Specifications</b>	<b>22</b>
<b>Warranty Statement</b>	<b>27</b>

# System overview

The ProCair Prime Advanced is an alternating mattress replacement system for the prevention and treatment of skin breakdown and pressure injuries in patients of high to very high risk. It is designed to replace your existing bed mattress on either a standard or profiling electric bed frame.

The revolutionary Universal Therapy System™ (UTS) features an integrated control unit within the mattress, cyclically inflating and deflating the air cells in an alternating pattern, to help prevent arterial and venous capillary occlusion in the patient's surface tissue - maintaining the flow of blood and lymphatic fluids through these tissues to provide essential oxygen and remove metabolic waste.

A full-size, movable foam sheet can be positioned at the bottom of the mattress, working in conjunction with the cell in cell air cells to provide bottoming-out protection. Alternatively, the foam sheet can also be positioned on top of the air cells to convert the mattress to a hybrid support surface, providing additional comfort for the user as well as gentle alternating therapy.


The ProCair Prime Advanced can be easily operated by the caregiver as well as the patient (subject to advice by a healthcare professional) via a simple-to-use handset with a long, flexible cord.



The system consists of the following components:

- ~ Mattress replacement with integrated control unit, control handset, power cord and CPR release
- ~ User Guide
- ~ Carry bag

It is recommended that all packing materials and the User Guide be kept in the carry bag, for ease of storage and/or transport.






# Safety precautions

The purpose of the following safety precautions are to direct attention to possible dangers. The safety symbols and their explanations require careful attention and understanding.




The safety warnings by themselves do not eliminate any danger. The instructions or warnings they give are not substitutes for proper accident prevention measures.



For your own safety and the safety of equipment, always take the following precautions.

## General Safety Precautions


-  Read all instructions before using this medical device
-  This system must be used on top of an appropriate sized bed frame and the appropriate operating environment as stated in this User Guide.
-  Before commencing set up or installation, ensure the power is switched off and disconnect the power cord from the control unit.

Novis Healthcare recommends using the cord retention loops on the mattress base where possible and attaching it to an electrical outlet by the head of the bed.




-  Minimise layers between patient and mattress and secure bed sheets loosely so as not to affect the alternating cell movement. As part of a sensible pressure injury prevention strategy, avoid wearing clothing that may cause areas of localised damage due to creases, seams, objects in pockets, etc.
-  Never use sharp objects or electrically heated blankets on or under the system.
-  Product top cover may present a suffocation risk. It is the responsibility of the caregiver to ensure that the patient can use this product safely.

-  Avoid blocking the air intakes of the control unit, located on the right side of mattress base (when viewed from foot end).
-  Bed frames used with the systems can vary greatly depending on the specific healthcare setting (ie. hospitals, aged care, home care, etc).

It is the responsibility of the caregiver to take the necessary precautions to ensure the safety of the patient. This includes, but is not limited to, the appropriate use of side rails to prevent falls.

-  Only the control unit and mattress combination as indicated by Novis Healthcare should be used, otherwise the correct function of the product cannot be guaranteed.

## User Capacity




-  The maximum recommended safe working load for this system is 200 kilograms (440 pounds).
-  The minimum recommended safe working load for this system is 30 kilograms (66 pounds).
-  Do not exceed this safe working load or you risk injury to the patient or carer and damage to the product.

# Safety precautions




## Protection Against Hazards

### Fluids


Avoid fluid contact on any part of the control handset or control unit, by keeping mattress top cover zipped closed and the waterfall skirt down. If spills do occur:

- ~ Turn off control unit power and disconnect the unit from mains electricity supply.
- ~ Immediately clean fluids from the casing by wiping with a soft cloth.
-  Ensure there is no moisture in or near the power inlet, control handset and power cord before reconnecting the power supply.
-  Check the operation of controls and other components around the spill area.
-  Fluid or liquid remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and carers.

### Explosion Hazard

-  Equipment is not suitable for use in the presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide.
-  Do not use in the presence of smoking materials or open flame – air flowing through the mattress will support combustion.
- Do not open the control unit – risk of  electrical shock. Refer servicing to qualified service personnel.




### Disposal

-  Dispose of all components (control unit including batteries, air filter, air cells, mattress cover and base) according to local procedures and regulations or contact Novis Healthcare for advice.

### Power Cord

The power cord and pump shall be placed at the foot side of the patient to prevent any risk of strangulation due to the cord.

Periodically inspect the power cord for damage.

-  The system should never be operated with a worn or damaged power cord.
-  Keep the cord away from heated surfaces and places where it may be difficult to unplug the device, such as shelves, wardrobes, or ceiling sockets.
-  Should the power cord be found to be worn or damaged, contact Novis Healthcare for a replacement.

### Interference

Although this equipment conforms to the intent of directive IEC 60601-1-2<sup>1</sup> in relation to Electromagnetic Compatibility, all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact Novis Healthcare.

### Health hazard

If pain, irritation, numbness, swelling, or redness occurs, discontinue use and contact a healthcare professional.

1 IEC 60601-1-2: Medical Electrical Equipment - Part 1: General Equipments for Safety, Amendment No. 2. Collateral Standard. Electromagnetic Compatibility Requirements and Test).

# System preparation

Carefully unpack the system and inspect each item for any damage that may have occurred during transit and handling. Any damage or missing components should be reported to Novis Healthcare as soon as possible.

Prior to usage, please allow the product to acclimate to ambient conditions for at least 30 minutes if the ambient temperature is 20°C. This helps ensure optimal performance and prevents potential damage caused by sudden temperature changes from the storage temperature range of 5°C to 60°C.

**⚠** Confirm there are no sharp objects in the immediate area which may risk damage to the mattress replacement.

1 Remove existing mattress from bed, release the main strap to unroll ProCair Prime Advanced mattress. Place the mattress on top of the bed base – printed top cover facing upwards.



2 Attach to the bed by securing the adjustable straps, located on the underside of the mattress base under each bed end.


On a profiling bed, secure the straps around the moveable sections of the base. Ensure the buckles are securely fastened and straps are pulled tight.

**⊗** Do not secure mattress straps to bed side rails – straps will tear.



**⚠** Ensure that straps do not interfere with the operation of the bed, and that the mattress is properly secured. Failure to do so could result in patient injury or equipment damage.

3 Check CPR sealing valve is closed – the turning tab and the arrows must be aligned to 'CLOSED' position.



# System preparation

- 4 Check the internal quick release control unit air connectors are securely connected.

Open the top cover by unzipping the right side of the mattress (zipper located at foot end), check the two connectors are secure by pushing the male and female ends together, there should be no movement.

If a connection is open, a click will be heard once connector is firmly closed.



- 5 Snap the hanging hook over the back of the control handset as shown. Hang the control handset over the foot end of the bed.

Ensure the handset is secure and not at risk of falling off before use to prevent equipment damage.




- 6 Feed power cord through the cord retention loops along either side of the mattress base. Insert power cord plug into the side of the control unit, then connect to an appropriate electrical outlet and switch on mains power.

The Standby indicator on the handset will glow amber, confirming the control unit is connected to a power source.



- ⚠** Ensure the power cord is not under strain; is free from obstruction; and is secured safely so as not to be a trip hazard.

# System preparation



7 On the control handset, press and hold the Power button for a minimum of three seconds.

The Power indicator will glow green to indicate the system is operational and automatically inflating.

While inflating, the mode/pressure indicator will alternate between "oo" and "00", indicating mattress is being inflated to maximum pressure in preparation for use. Allow up to 45 minutes for complete inflation.

**⚠** Do not lie a person (or any weight) on the mattress during initial inflation.

8 When initial inflation is complete, a beep will be emitted and the mode/pressure indicator will display "AL" for **AL**ternating, to indicate the system is ready for use in Alternating Mode.

**⚠** The system defaults to pressure setting 03 upon full initial inflation. Press "AUTO SET" to activate automatic weight detection. This feature monitors patient weight every alternation cycle and automatically selects the appropriate pressure setting accordingly.

This feature can be deactivated for manual pressure setting (see page 14 for details).



Dynamic



Hybrid

## Hybrid / Dynamic configuration

9 The ProCair Prime Advanced is **initially configured as a full dynamic mattress**, with the patient in direct contact with the air cells for maximum therapeutic benefits.

To convert to hybrid mode, refer to instructions on page 13.

# System preparation

## Universal Therapy System™ Mattress Configuration



### Dynamic Mode

The ProCair Prime Advanced is initially configured as a full dynamic mattress, the patient is placed in direct contact with the air cells for optimum pressure offloading during alternation, and the foam sheet is positioned in its storage compartment at the bottom of the mattress to provide bottoming out protection.

### Hybrid Mode

Hybrid mode allows greater comfort for the patient by adding a softly padded foam sheet over the air cells to allow immersion of the patient's body over the mattress surface, while providing a limited degree of alternation and pressure offloading (in comparison to dynamic mode).

#### To convert to hybrid mode

- 1 Unzip zippers from the foot end of the mattress and remove the top cover.
- 2 Undo 4 press studs on one side of the mattress base cover
- 3 Remove the foam sheet from its compartment at the bottom of the mattress.
- 4 Align the 7 press studs on each side of the foam sheet with those along the top edges of the mattress base, fasten the foam sheet over the air cells.
- 5 Re-zip the top cover in place



#### To return to Dynamic Mode

Simply remove foam sheet from the top of the mattress by undoing the press studs and returning the sheet to its storage compartment at the bottom of the mattress. If necessary, the foam sheet may be cleaned and dried before returning to its compartment, the foam is NOT washable.

- The patient must be removed from the mattress during the conversion between Hybrid and Dynamic modes.
- Consult your health professional for the suitability of Hybrid Mode according to the patient's condition and pressure injury risk.

# Intended use

## Indications

The ProCair Prime Advanced mattress replacement systems are indicated for:

- ~ The prevention and treatment of skin breakdown and pressure injuries in patients of high to very high risk.

## Contraindications

Patient conditions for which the application of pressure therapy on the ProCair Prime Advanced mattress replacement systems contraindications include:

- ~ Instable spinal cord injury
- ~ Cervical traction

## Intended care setting

Intended care settings for the ProCair Prime Advanced mattress replacement systems are:

- ~ Home healthcare
- ~ Professional healthcare

## Working environment

- ~ Temperature: 15°C to 35°C (59° F to 95° F)
- ~ Humidity: 30% to 75% non-condensing

## Shipping / storage environment

- ~ Temperature: 5°C to 60°C (41° F to 140° F)
- ~ Humidity: 15% to 90% non-condensing

## Connecting system to other devices

The ProCair Prime Advanced single mattress replacement can be fitted to most standard hospital or single bed bases.

The ProCair Prime Advanced king single mattress replacement can be fitted to most king single hospital or king single bed bases.

Therapeutic devices should only be used in accordance with manufacturer's instructions and under the consent, supervision and management of a suitably qualified health professional.

Novis Healthcare accepts no liability for any use, change or assembly of the product other than that stated in this User Guide. Refer to our Warranty Statement for more details.

## Intended user profile

The user (operator) of this device must have at least eight years of education and have the ability to read and understand English and Westernized Arabic numerals.

# Patient set-up

Once initial inflation is complete, a patient may be placed on the system.


1 Once the mattress is fully inflated, bedding can be replaced.

Secure sheets loosely enough to ensure they do not interfere with cell alternation.

2 Place the patient on to the ProCair Prime Advanced mattress.

The system defaults to pressure setting 03 upon full initial inflation. Press "AUTO SET" to activate automatic pressure setting. This feature constantly monitors patient weight (30-200kg / 66-440 lbs) and automatically selects the appropriate pressure setting accordingly in real time.

3 Perform a 'bottoming out' test (a test to ensure the patient is adequately suspended away from the base).



## 'Bottoming Out' test

1 Check system is in alternation mode by ensuring the mode/pressure indicator is displaying "AL" and that one set of air cells is inflated while the other set is deflated.

You may need to unzip the cover to feel the cells for inflation.

2 With the patient lying supine, unzip one side of the top cover just past sacral region (lower spine).

3 Slide your hand underneath the patient and feel for a deflated cell under the patient's lower spine (in Hybrid Mode, place hand underneath the foam sheet under the patient's lower spine) your hand should easily slide between patient and base.


4 If your hand can pass under the patient, then patient is adequately suspended. If not, manually adjust pressure to 'firmer' and wait at least one cycle (14minutes) for pressure to increase before repeating step 3.

If manual pressure adjustment fails, press Max Inflate to force mattress to full inflation. Wait at least one cycle (14 minutes) for pressure to reach maximum pressure, then press Alt to return to an alternation cycle. Wait at least one more cycle (14 minutes) for pressure to increase before repeating step 3.

5 We recommend repeating the Bottoming Out test at least 14 minutes after any manual pressure readjustment.

# Operation

## Control handset



### B Power

Press and hold the Power button for at least 3 seconds to turn the system power on and off.

Power LED (**A**) - in operation

Standby LED - power connected

### C Battery in use

Battery in use LED indicates the device is using battery power when illuminated in green.

### D Battery status

<b>Green LED lit</b>	The battery is fully charged
<b>Green LED flashing</b>	The battery is charging
<b>Amber LED flashing</b>	The battery power is low
<b>LED is off</b>	The battery is malfunctioning

### F Mode/pressure indicator

Normally displays the current mode of operation (AF, St, AL, FL). Also displays pressure levels from 01 to 05 for 10 seconds when mattress pressure is manually adjusted, before returning to Mode display.

### G Manual pressure select

After initial inflation, the system defaults to setting 03 (40mmHg), manual adjustments may be made according to patient weight, levels of support and comfort required.

Press button to cycle pressure settings from 01 to 05, pressure levels are as follows:

<b>01</b>	20mmHg
<b>02</b>	30mmHg
<b>03</b>	40mmHg
<b>04</b>	50mmHg
<b>05</b>	60mmHg

### H Auto pressure setting

Press to activate and deactivate automatic pressure setting. Once activated, as indicated by the green LED above (**E**) the system will automatically detect patient weight and constantly adjust pressure accordingly in real time as necessary.

## Modes of operation

### I Mode selection

Press button to cycle between modes, as indicated by the mode/pressure indicator:

AF	Auto Firm (Max Inflate)
St	Static
AL	Alternating
FL	Fowler Mode

Refer to the next page on details of the above modes of operation.

### K Alert indicator

LP	Low pressure due to possible air leak
PF	Power failure, mains power disconnected
Sd	Sensor detect

For further details refer to **Troubleshooting** on Page 20

### L Alert mute/panel Lock button

If alert sounds, press once to mute alert buzzer for 20 minutes to minimise disturbance.

**Panel Lock** - Long press for 3 seconds to lock control panel to prevent tampering or unintended operation, an amber LED indicator above will illuminate (**J**)

To unlock panel, long press for 3 seconds again. After 5 minutes of inactivity, the panel will re-lock automatically.

In **Alternating Mode**, alternate mattress cells inflate and deflate following a fixed cycle time of 14 minutes, with the exception of static head cells.



Alternating mode is used for normal therapeutic function.

In **Static Mode**, all mattress cells are inflated thereby suspending the alternation cycle and providing constant low pressure therapeutic patient benefits, also to facilitate nursing or feeding.





In **Auto Firm Mode**, all mattress cells are inflated to maximum pressure to create the firmest possible base for stable patient handling and transport or other special circumstances.




In Auto **Fowler Mode**, the mattress pressure is increased to accommodate extra load in the sacral region. Auto Fowler Mode engages when the bed backrest is raised and the mattress is inclined to 30° or more. This feature helps to prevent the risk of bottoming out.



 The system will operate in Auto Firm Mode for a maximum of 20 minutes, after which it will automatically revert to Alternating Mode for patient safety. The system will operate in Static Mode until manually changed to another mode.

 During power failure/outage, the ProCair Prime Advanced Universal Therapy System continues to provide dynamic alternation therapy for up to 5 hours with automatic battery back up. It may take up to 16 to 20 hours to fully charge a depleted battery. **PF** alert code will be displayed on the Alert Indicator, accompanied by an audible alert. The system will return to its normal operation when mains power supply is resumed.

 When using a pump with the battery in a normal situation, if a pump is turned off immediately after losing the AC power, it might be the battery is out of power or failed. In this case, connect the pump to the AC power to charge the battery for 16-20hours. If the pump still cannot be powered on by the battery, it means that the battery is out of order. Please contact a qualified and authorised professional to replace the battery.

# Operation

## Quick Twist CPR

Rapid deflation of the mattress may be required for emergency treatment (or to pack up the unit).

The Quick Twist CPR valve is located at the top of the mattress, to the right of the patient's head.

If emergency resuscitation is required, turn the CPR valve so the arrows are aligned with the 'OPEN' position markers. This will rapidly deflate the entire system, including static head cells.



CPR in the closed position.  
Twist in either direction to open

To reinflate the system after the Quick Twist CPR valve has been released, turn the tab and align the arrows with the 'CLOSED' position markers, ensure control unit is switched on and wait for the system to regain optimal pressure.



CPR in the open position.  
Twist in either direction to close

## Transport function



To prepare for patient transport, press the Mode button and select Static mode, wait 14 minutes to ensure all cells are fully inflated.

Ensure Quick Twist CPR valve is set to the CLOSED position. Switch off control unit by long pressing the Power button for 3 seconds.

Open top cover, disconnect the air connectors from the control unit and connect the male connector to the female connector to ensure air does not escape from the mattress. Tuck hoses under an air cell.



Air will remain in the system for up to 24 hours, depending on patient and environmental circumstances.

## Deflation and storage

- 1 Press the power button for a minimum of **three seconds** to switch off the control unit.
- 2 Switch off mains power and unplug the power cord from the mains outlet.
- 3 Turn the Quick Twist CPR to OPEN to release air and deflate all cells.
- 4 Once air has been released from the system, detach the mattress from the bed by unfastening all mattress straps, place handset and power cord on the centre of the mattress, then fold and roll the mattress from head end to foot end.

- 5 Fasten main strap to keep mattress in rolled position. Return all items to the custom carry bag for safe keeping.

- 6 The bag is designed to fit snugly around the mattress. It may be easier to turn the open carry bag upside down, and then wrap the bag around the rolled up mattress, turn over bag, then tuck in the mattress as you zip up the bag.






Power button




Main strap



# Care and cleaning

-  To prevent cross-contamination, the mattress should be examined and disinfected between patient use, and a liner or bed sheet is recommended for added protection.
-  Clean the mattress in accordance with local infection control policy and government regulations. Failure to do so could cause patient or personal injury.
-  The mattress is not protected against excessive amounts of fluid. Do not immerse the control unit in fluid.
-  Switch off and disconnect the control unit from mains power supply before cleaning. Failure to do so could result in equipment damage or electric shock.
-  Do not use high temperature autoclave steam cleaning devices or phenolic based products for cleaning. This could result in damage to the equipment and may result in damage to the polyurethane coating, or negate the biocompatibility properties of the fabric.



## Cleaning and infection control

-  It is recommended that the ProCair Prime Advanced system is cleaned **every two weeks** if in constant use.

### Top cover cleaning

Unzip and remove the top cover from the base before washing (refer page 19 for instructions).


For basic care and cleaning, wipe down with warm water containing PH neutral detergent. The top cover can also be machine washed at a maximum of 95° C (203° F) using neutral detergents.

-  Refer to the top cover wash tag for detailed cleaning instructions.
-  Do not use system without top cover.

### Base and foam sheet cleaning

Swab the mattress base, air cells and the foam sheet with a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine). Dry thoroughly before reassembly.

Do not machine wash or tumble dry the air cells, mattress base or the foam sheet.



-  If cleaning or disinfection is required, do not allow fluid to enter air cells and air hoses.

### Control Unit /Handset cleaning

Disconnect control unit from mains power before cleaning. Gently wipe down the external case with a soft cloth.

Soak the cloth in warm water containing mild PH neutral detergent, and wring any excess water before gently wiping all external controls. Repeat the process with a dry cloth to remove excess moisture. A soft bristled nylon brush can be used to gently clean crevices.

Ensure the control unit is disconnected from mains power before cleaning.

-  Do not spray disinfectant directly on to the control unit, or immerse the unit in water or
-  other fluid.

### Disinfection

The mattress, top cover and control unit may be decontaminated by a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine). Dry thoroughly before use.

For infection control, swab with a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine). Dry thoroughly before reattaching and use.

### Top Cover Removal

- 1 Raise the waterfall skirt and locate the zippers at the foot end of the mattress.



- 2 Starting with either zipper, run the zipper along the side of the mattress towards the centre of the head end.



- 3 Repeat with other zipper. The top cover can now be detached from the mattress base.


To reattach the top cover to the mattress base, reattach the zippers, then run them along the sides of the mattress back towards the centre of the foot end.



# Troubleshooting


## Alert Codes

An audible alert, accompanied by an Alert Code display, indicate the control unit or mattress is experiencing a fault. The Alert Code display will continue to flash until the fault is cleared.

The audible alert can be silenced for 20 minutes by pressing the Alert Mute button. It will reactivate if the fault is not rectified or if a new fault is detected.




Alert Code	Trigger	Solution
 <b>LP = Low Pressure</b>	Air cells have failed to reach the pre-set pressure after two cycles	<p>Ensure mains power is turned on and power cord is connected to mains and control unit.</p> <p>Check control unit/mattress air connections are fitted securely.</p> <p>Ensure control unit is turned on.</p> <p>Ensure CPR valve is set to CLOSE position. Replace CPR valve if air leak is found.</p> <p>Check air intake from filter is not blocked by linen/dust. Replace with new filter if needed.</p>
 <b>PF = Power Failure</b>	Power Failure	<p>Check control unit is connected to the mains power supply.</p> <p>Check for loose connection on plug and main power is switched on.</p> <p>Check condition of power cord and plug. Check if mains socket is faulty.</p>
<b>Control unit does not operate; no display or lights</b>		<p>Check control unit is connected to the mains power supply.</p> <p>Check for loose power cord connection and ensure main power is switched on.</p> <p>Check condition of power cord and plug. Check if mains socket is faulty.</p>

 For faster mattress reinflation once the air leak has been closed, select Auto Firm Mode (**AF**) and wait until the LP alert code stops displaying.

Select Alternation Mode (**AL**) to resume alternation.


Alert Code	Trigger	Solution
<b>Patient is sinking or 'bottoming out' whilst lying flat on the mattress</b>		<p>The pressure may be set too low for the patient's weight – increase the pressure setting by manually setting a higher pressure level by pressing the "MANUAL" button.</p> <p>Check for air leaks in the mattress and air hoses.</p>
<b>Control unit controls lock up, or freezes</b>		<p>Turn off and unplug the control unit from mains power.</p> <p>Rest the control unit for one minute before reconnecting the control unit to mains power and switching on.</p>
<b>Air flow output varies unstably or stops erratically</b>	<p>This may be caused by electro-magnetic disturbance or unstable mains power supply.</p>	<p>Ensure system is used in an area with a stable power supply, if this is not possible, consider connecting an Uninterruptible Power Supply (UPS) unit or a battery power/backup between the system and the mains power supply to stabilise power.</p> <p>Please ensure the battery has sufficient power while using battery power/backup.</p> <p>Restart system by powering OFF, and then powering ON again.</p>

 If the problem persists, move patient to another mattress and contact Novis Healthcare.

# Technical Specifications

Prime Advanced		Single	King Single	
System	System code	APMPC-R03A	APMPC-R03AK	
	Capacity	30 - 200 kg		
	No of cells	21, including 3 static head cells and 4 narrow foot cells		
	Compliance	IEC60601-1, IEC60601-1-2 and IEC60601-1-11		
	ARTG	289458		
Control unit	Part number	APMPC-CU07		
	Control system	Digital micro controller		
	Casing material	ABS		
	Cycle time	14 minutes (fixed)		
	Dimensions	L 295 x W 140 x H 95 mm		
	Voltage /current	AC100 - 240V / 50Hz-60Hz / 0.2-0.1A		
	Fuse Rating	T2AH 250V		
	Power Cable	5 meters, non-shielding		
	Operation Range	Atmospheric Pressure 700 to 1060 hPa		
	Min / max pressure	20 to 60 mmHg +/- 6 mmHg		
	Protection type	Class II Type BF (Not AP or APG Type)		
	Ingress protection rating	IP22	SOLIDS: Protected from being touched by fingers and objects greater than 12.5 mm. LIQUIDS: Against water: Dripping water when tilted at 15°.	
	Handset casing material	ABS		
	Handset cable	3.5 meters (straightened)		
	Handset part number	FO-NVS1006		
Mattress dimensions	Length	2000 mm	2000 mm	
	Width	900 mm	1050 mm	
	Height	200 mm	200 mm	
	Weight	14 kg	15 kg	
	Material	Top cover	PU laminated nylon	
		Base cover	PVC laminated polyester	
		Air cell	TPU laminated nylon	
Operating environment	Air humidity	Operation	30% to 75% non-condensing	
		Storage	15 % to 90% non-condensing	
	Ambient temperature	Operation	15° C to 35° C (59° to 95° F)	
		Storage	5° C to 60° C (41° to 140° F)	

## Waste disposal



This product has been supplied from an environmentally aware manufacturer that complies with the European Community's Waste Electrical and Electronic Equipment Directive (WEEE).

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to the legislation. Please be environmentally responsible and contact your local authority on available options to recycle this product at its end of life.

## Service life

The service life of the internal electrical power source (battery) is 2 years.

The expected service life of a control unit and a mattress is highly dependent on frequency of use, servicing, care and maintenance.

To maintain the condition of the alternating mattress system, service the system regularly according to the schedule recommended by Novis.


## Maintenance


The device should only be repaired, including the installation or replacement of the battery pack, by a qualified and authorized professional. The product has no user/operator serviceable parts.


Preventative maintenance needs to be conducted every 12 months. Please contact Novis at **1300 738 885** or a certified maintenance technician.

Before performing maintenance or service, turn off the power, ensure the mattress is empty, and refrain from using the device.


Failure to do so may lead to exposure to hazards.


 Do not use unapproved accessories or attempt to modify, disassemble or otherwise misuse the ProCair Prime Advanced system or any of its components.

 All product specifications are subject to change without notice.

 Medical electrical equipment needs special precautions regarding EMC and needs to be installed according to the EMC information provided. Careful consideration of this information is essential when stacking or collocating equipment and when routing cables and accessories.

 RF mobile communications equipment can effect medical electrical equipment.

 X-rays, electric cauterizers, high-frequency or high-energy devices may cause interference with the device.

 Clip the core onto the power cable of the pump, approximately 70 cm from the pump, to prevent Radiated Emission (RE) from interfering with other medical electrical equipment.

# Technical specifications


The Procair Prime Universal Therapy System M31 Control Unit is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the Procair Prime Universal Therapy System M31 Control Unit should assure that it is used in such an environment.

## Guidance and manufacturer's declaration – electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PROCAIR PRIME UNIVERSAL THERAPY SYSTEM M31 CONTROL UNIT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliance	The PROCAIR PRIME UNIVERSAL THERAPY SYSTEM M31 CONTROL UNIT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

## Guidance and manufacturer's declaration – electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment –guidance (for home and professional healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ±8 kV	Contact: ±8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
	Air: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Air: ±2 kV, ±4 kV, ±8 kV, ±15 kV	
Electrical fast transient/ burst IEC 61000-4-4	± 2kV for power supply lines	± 2kV for power supply lines	Mains power quality should be that of a typical home and professional healthcare environment.
	± 1kV for input/output lines	± 1kV for input/output lines	
Surge IEC 61000-4-5	± 0.5kV, ± 1kV line(s) to line(s)	± 0.5kV, ± 1kV line(s) to line(s)	Mains power quality should be that of a typical home and professional healthcare environment.
	± 0.5kV, ± 1kV, ± 2kV line(s) to earth	Not applicable	
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0 % UT; 0,5 cycle 0 % UT; 1 cycle 70 % UT; 25/30 cycles	Voltage dips: 0 % UT; 0,5 cycle 0 % UT; 1 cycle 70 % UT; 25 cycles	Mains power quality should be that of a typical home and professional healthcare environment. If the user of the Procair Prime Universal Therapy System M31 Control Unit requires continued operation during power mains interruptions, it is recommended that the Procair Prime Universal Therapy System M31 Control Unit be powered from an uninterruptible power supply or a battery.
	Voltage interruptions: 0 % UT; 250/300 cycle	Voltage interruptions: 0 % UT; 250 cycle	
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The Procair Prime Universal Therapy System M31 Control Unit power frequency magnetic fields should be at levels characteristic of a typical location in a typical home and professional healthcare environment.
	50 Hz or 60 Hz	50 Hz and 60 Hz	

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance (for home and professional healthcare environment)
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz	3 Vrms: 0,15 MHz – 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the ProCair Prime Universal Therapy System M31 Control Unit including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts ( <b>W</b> ) according to the transmitter manufacturer and d is the recommended separation distance in metres ( <b>m</b> ).
	6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz	6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz	
Radiated RF IEC 61000-4-3	80 % AM at 1 kHz	80 % AM at 1 kHz	Interference may occur in the vicinity of equipment marked with the following symbol: 
	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	

**NOTE** UT is the a.c. mains voltage prior to application of the test level.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## Proximity to other electronic equipment

Recommended separation distances between portable and mobile RF communications equipment and the ProCair Prime Advanced Universal Therapy System M31 control unit

The ProCair Prime Advanced Universal Therapy System M31 control unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ProCair Prime Advanced Universal Therapy System M31 control unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ProCair Prime Advanced Universal Therapy System M31 control unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter ( <b>W</b> )	Separation distance according to frequency of transmitter ( <b>m</b> )		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in metres (**m**) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (**W**) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

# Technical specifications

## Manufacturer's declaration – electromagnetic immunity

### Test specifications for Enclosure Port Immunity to RF wireless communications equipment

Test frequency (MHz)	Band <sup>A</sup> (MHz)	Service <sup>A</sup>	Modulation <sup>B</sup>	Maximum power (W)	Distance (m)	Immunity test level (V/m)	Compliance level (V/m) (for home and professional healthcare)
385	380 – 390	TETRA 400	Pulse modulation <sup>B</sup> 18 Hz	1.8	0.3	27	27
450	430 – 470	GMRS 460, FRS 460	FM C ±5 kHz deviation 1 kHz sine	2	0.3	28	28
710	704 – 787	LTE Band 13, 17	Pulse modulation <sup>B</sup> 217 Hz	0.2	0.3	9	9
745							
780							
810							
870	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>B</sup> 18 Hz	2	0.3	28	28
930							
1,720							
1,845	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation <sup>B</sup> 217 Hz	2	0.3	28	28
1,970							
2,450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>B</sup> 217 Hz	2	0.3	28	28
5,240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation <sup>B</sup> 217 Hz	0.2	0.3	9	9
5,500							
5,785							

**NOTE** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

**A** For some services, only the uplink frequencies are included.

**B** The carrier shall be modulated using a 50 % duty cycle square wave signal. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

### Test specifications for Enclosure Port Immunity to proximity magnetic fields

The Procair Prime Universal Therapy System M31 Control Unit is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the Procair Prime Universal Therapy System M31 Control Unit should assure that it is used in such an environment.

Frequencies	Test Level (A/m)	Modulation	Dwell time (s)	Compliance LEVEL (V/m) (for home and professional healthcare)
30 kHz <sup>A</sup>	8	CW	3	8
134.2 kHz	65	Pulse modulation <sup>B</sup> 2.1 kHz	3	65 <sup>C</sup>
13.56 MHz	7.5	Pulse modulation <sup>B</sup> 50 kHz	3	7.5 <sup>C</sup>

**NOTE A** This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME AND PROFESSIONAL HEALTHCARE ENVIRONMENT.

**B** The carrier shall be modulated using a 50 % duty cycle square wave signal.

**C** r.m.s., before modulation is applied.

# Warranty statement

## This warranty is provided by

### Novis Healthcare

(ABN 45 102 735 491)  
of **Unit 12, 12 Mars Road Lane  
Cove West, NSW 2066.**

Novis Healthcare (Novis) products are manufactured to the highest quality standards and are thoroughly tested and inspected before leaving our factory. In addition to any statutory rights and remedies you may have, Novis warrants all of its products sold directly or via an Authorised Novis Australia Dealer against defective workmanship and faulty materials from the date of purchase by the end user for a period of thirty-six months unless otherwise specified for that product and its components.

## Warranty Claims

To claim under this warranty, please contact Novis Healthcare and have your receipt or proof of purchase available. Novis Healthcare may need to assess the defect before determining any claim,

and additional information may be requested to process your claim. Claims without proof of purchase may not be able to be processed.

Novis Healthcare may at its option inspect the goods on site or require them to be returned to its premises or one of its Authorised Service Agents in person or freight prepaid by you.

Novis will undertake at its option, to repair or replace, free of charge, each product or part thereof on the condition that:

- ~ The product found on examination, to be suffering from a manufacturing defect;
- ~ The product or relevant part has been serviced regularly by Novis or one of its Authorised Service Agents and has not been subjected to misuse, neglect or been involved in an accident;
- ~ The repairs are not required as part of normal wear and tear.

At our option

- ~ Goods repaired may be replaced by refurbished good of the same type rather than being repaired.
- ~ Refurbished parts may be used to repair goods.

Novis Healthcare will not be held responsible for any repair other than those carried out by it or one of its Authorised Service Agents.

Warranty repairs do not extend the length of the warranty period.

## Limited Liabilities

Our liability under this manufacturer's warranty is subject to us being satisfied that a defect was caused by faulty parts, manufacture or workmanship, and was not caused or substantially contributed to by other factors or circumstances beyond our control, including (but not limited to) defective installation, maintenance or repair, product modification or alteration, any neglect, misuse, or excessive use, normal wear and tear or failure to follow manufacturer's instructions.

## Important notice for Australian consumers:

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law.

You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. To obtain compensation, you will need to provide documentary evidence of the loss or damage suffered and documentary evidence that such loss or damage was a reasonable foreseeable consequence of a failure Novis Healthcare to comply with a consumer guarantee under the Australian Consumer Law.

Subject to the provisions of the Australian Consumer Law, Novis Healthcare excludes, to the fullest extent permitted by law, all liability in respect of loss of profit or other economic loss, direct to indirect or consequential, special, general or other damages or other expenses or costs which may include negligence.







1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au



1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service



1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au