

Proper battery care is essential to ensure your hoist remains safe, reliable, and ready for use when needed. Following these simple guidelines can help maximise battery life and reduce the risk of unexpected equipment downtime.

DAILY USE

- Always return the hoist battery to its charger when not in use.
- Ensure the charger is plugged in and switched on at the power outlet.
- Check that the charging indicator light is illuminated.
- Avoid allowing the battery to become completely flat before recharging.

CHARGING GUIDELINES

- Charge the battery in a clean, dry environment after every use.
- Allow the battery to charge fully whenever possible.
- Do not disconnect the charger before charging is complete unless the hoist is required for use.
- If the hoist is not being used regularly, leave it connected to the charger to maintain battery health.

BATTERY SAFETY

- Only use approved chargers supplied for your hoist.
- Do not attempt to open, modify, or repair the battery.
- Keep batteries away from water, excessive heat, and naked flames.

SIGNS YOUR BATTERY MAY NEED ATTENTION

Contact your service provider if you notice:

- Reduced lifting capacity or operating time.
- The hoist requires frequent charging.
- Charging indicator lights do not operate correctly.
- Battery warning alarms activate regularly.
- Physical damage to the battery or charger.

NEED ASSISTANCE?

If you have any questions regarding your hoist, battery care, servicing, or repairs, please contact:

Rehab Hire & Sales | Repairs and Service

†: 1300 000 030 and choose option 4.