



**HEAD LIFT & KNEE BREAK**

Head Lift 62°, Knee Break 42°



**200KG SAFE WORKING LOAD**

Queen & King 200kg per side



**5 YEAR WARRANTY**

Refer to our Warranty Card for details



**BASE HEIGHT**

37cm from floor to top of base  
Non Height Adjustable

The IC111 Bed is our low risk extra mobility option. These domestic looking bases offer head lift and knee break. Queen and king size are available with split mattresses with no gap and individual controls.


Choice of Mattresses

Available in Onyx & Stone fabrics

Headboards - optional extra



Long Single  
Base:  
(2150 x 960)  
Mattress:  
(2030 x 900)



King Single  
Base:  
(2150 x 1120)  
Mattress:  
(2030 x 1070)



Queen (per side)  
Base:  
(2150 x 790)  
Mattress:  
(2030 x 750)



King (per side)  
Base:  
(2150 x 960)  
Mattress:  
(2030 x 900)





## IC111 Headboards

Upholstered headboard to match IC111  
Available in all bed sizes  
Easily attachable and detachable



## Companion Base

Non Adjustable

Non mechanical upholstered slat base  
Available in onyx and stone fabric  
To sit beside your homecare bed

Long Single  
Base: 2150 x 960  
Mattress: 2030 x 900

King Single  
Base: 2150 x 1120  
Mattress: 2030 x 1070

## Accessories



### BEDSTICK

Removable from brackets  
Height: 50cm from top of base



### HIGH SIDE RAIL

Removable from brackets  
Height: 72cm from top of base  
Width: 36cm



### LOW SIDE RAIL

Removable from brackets  
Height: 42cm from top of base  
Width: 68cm



### FULL LENGTH FOLD-DOWN RAIL

Fold down mechanism  
Height: 34cm from top of base  
Width: 138cm  
Available in left and right hand side



### OVER BED POLE

Free standing  
Overall Height - 175cm  
Adjustable handle height