

Elite Deep Back



Outer cover:

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

Back shell:

- Stylish, lightweight and breathable aluminium back shell with 6" contour depth

Hardware:

- "Easy Set" hardware installed and adjusted in minutes - adjustments can be made while user is seated in the wheelchair

Accessories

Mounting Pins



Chest Harness Bracket



Privacy Flap



Split



Heavy-Duty Easy Set



Extended Range



Fixed Space Saver



Tech Specs

Elite Deep Back	 Height 12", 16", 18", 20"	 Width 14" - 20"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40° range***	 Product weight from 1.3kg (3lb)	 Weight Capacity 136kg (300lb)	 Limited Warranty 2 years
Elite Deep Back Heavy Duty	Height 16, 20"	Width 21 - 24"	Depth adjustment 2" forwards, 2" backwards	Angle adjustment 40° range***	Product weight From 2.7kg (6lb)	Weight Capacity 227kg (500lb)	Limited Warranty 2 years

* 14" height only available up to 18" wide

*** Limitation of depth/angle adjustment if used on chairs narrower than model width



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®