



THE KNEE SCOOTER

**KEEPING YOU
MOBILE**



KC2

KNEE SCOOTER

KC2 Knee Scooter features a comfortable cushioned knee support, front wheel rotation, and safety hand brake which can also be locked for stability. Its easy to fold design makes KC2 the ideal choice for transport and storage.



TECHNICAL SPECIFICATIONS

Overall Length	77cm
Overall Width	42cm
Overall Height	Up to 114cm
Frame Material	Mild Steel
Weight	10.8kg
Handlebar Adjustment	90cm - 114cm
Knee Support Adjustment	40cm - 57.5cm
Braking	Handbrake Standard
Front Wheels	8" Solid Polyurethane
Rear Wheels	8" Solid Polyurethane
Weight Capacity	158kg
Warranty	1 year components / Lifetime on the frame

KC2 AT KNEE SCOOTER

KC2 AT, All Terrain Knee Scooter is a popular choice for comfortable, stable, portable and durable non weight-bearing rehabilitation support. Made of durable steel construction with 10' Pnuematic tyres and, the KC2 AT Knee Scooter is robust. The height-adjustable handle post and ergonomic handlebar design make the KC2 AT Knee Scooter a comfortable ride.



TECHNICAL SPECIFICATIONS

Overall Length	75cm
Overall Width	42cm
Overall Height	Up to 115cm
Frame Material	Mild Steel
Weight	10.5kg
Handlebar Adjustment	92cm - 115cm
Knee Support Adjustment	10" Pneumatic Tyres
Braking	Disk Brakes
Front Wheels	10" Pneumatic Tyres
Rear Wheels	10" Pneumatic Tyres
Weight Capacity	158kg
Warranty	1 year components / Lifetime on the frame

ACCESSORIES

SHEEPSKIN KNEE PAD COVER



- Improves knee pad comfort.
- Reduces the body pressure and friction at points of contact.
- Improves breathability. Cooler on hot days and warmer on cold days.
- Elasticised real sheepskin cover easily installs over any knee scooter model. Material lined for additional comfort.

BOTTLE HOLDER

Multi-purpose cup or bottle holder suitable for Knee Scooter, Walker, Stroller, Wheelchair or bike.



It is lined inside with silver insulated to keep your drinks cold and a riveted hole to drain any spilt liquids – Mesh pocket in front is suitable for a phone or wallet.

Height	10.5kg
Material	Black PE Material
Diameter	10cm

A COMPACT MOBILITY
AID TO ASSIST PEOPLE
WHO HAVE A LOWER
LEG INJURY KEEP ACTIVE.
KEEPS YOUR INJURED
FOOT NON-WEIGHT
BEARING WHILE YOU
MOVE ABOUT FREELY.



SAFETY

Replaces crutches, which can be unstable, difficult to use and can increase the chance of re-injury. The Knee Scooter is a stable base, easy and enjoyable to use. Suitable for all ages.

COMFORT

Fully adjustable for comfort. No more sore armpits and arms caused by using crutches.

INDEPENDENCE

Lightweight and compact, the Knee Scooter makes it easy to manoeuvre around the home, office, restaurants or shopping with friends. The Knee Scooter can get you there with ease, even in tight spaces.

FASTER RECOVERY

Improves mobility and optimises independence in the community and at work.



 knee scooters
A U S T R A L I A

Phone 0432 847 576

kneescooters.com.au

