## Start M2 - Start 4 configurable model

The Start M2 is a fully configurable, high-quality wheelchair model for which all of the options of the Start range are available. This model offers the maximum number of combination possibilities with regard to dimensions, colours and adjustments which makes it possible to provide even large and very large users with an optimal device. Drum brakes are also available options, as are elevating footrests or amputation legrests. The Start M2 also has an adaptive component. The rear wheel adapter allows you to adjust the wheelchair precisely to the user's needs.

There are many options available with the new Start M2 for the fitting of stroke patients with hemiplegia and cerebral palsy patients. The Start M2 is just the right choice for custom fittings in the lightweight segment.



## Technical data

User Weight	130 kg
Seat width	380 / 405 / 430 / 455 / 480 / 505 mm
Seat depth frame size 1	405 / 430 / 455 mm
Seat depth frame size 2	455 / 480 / 505 mm
Back height	400 / 425 / 450 / 475 mm
Anterior seat height	380 – 510 mm
Posterior seat height	380 – 510 mm
Seat angle	selectable
Lower leg length	250 – 470 mm (continuously adjustable)

## Product highlights

- Fully configurable lightweight wheelchair
- Patented, tool-free, quick seat depth adjustment (up to 50 mm)
- Very low weight with full Ottobock quality (from 14.7 kg total weight, from 9.6 kg transport weight\*)
- Antibacterial push handles\*\*
- Excellent component compatibility
- New adjustable economy back upholstery
- Larger centre of gravity adjustment range
- Larger selection of seat depths, including for people with a height over 2 m
- Larger selection of seat widths
- More options for height adjustment and seat tilt
- · Wider foot space for hemiplegics, which allows "foot propelling"
- New one-hand drive for hemiplegics
- Numerous updated and improved accessories
- · Two differnt frame lengths for optimal seat depth adjustment
- Selection of 6 different colours (anthracite metallic, jet black RAL 9005, ice blue, silver metallic, saturn blue metallic, white RAL 9016)
- Side panel with armrest, height can be adjusted (without tools) or set (with tools), adjustable depth

\* The transport weight is the total weight of the wheelchair minus the weight of the rear wheels, footrests and seat cushion.

\*\* The antibacterial effect of the push handles has been verified by an independent institute.