

## THE MOST STABLE CUSHION FOR WHEELCHAIR USERS

### Trochanter support

firmer density to provide lateral stability | redistributes pressure to help create an ischial well

### Ischial well - Two compartments

for superior pressure reduction in the IT-coccyx area & extra positioning possibilities

Mid height zip closure for easy adjustment when client is seated | cushion keeps its shape, enhancing the user's skin protection

Front compartments provide leg alignment and stability | easily adjustable for individual positioning needs

**PRESSURE INJURY TREATMENT PLAN**

**100% BREATHABLE**

**SOLID SHAPE**

**MACHINE WASHABLE**

**READY FOR USE**

**145 STANDARD SIZES**

- ▲ Nine adjustable compartments for optimal stability
- ▲ Expanding the user's reach
- ▲ Extra positioning of the thighs

The Vicair Vector 02 is a pre-contoured cushion which consists of 9 individual compartments. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin at the ischials and coccyx by re-distributing pressures, stabilizing the pelvis and positioning of the femurs. For specific needs the filling of each of the 9 compartments can be easily adjusted. Allowing you to change the density of each compartment benefitting users who: (1) Need correction or accommodation of an asymmetric pelvis. (2) Want (partial) off-loading of the ischial and coccyx area. (3) Need correction or accommodation for better thigh alignment. (e.g. abduction, one side foot propulsion, opening pelvic thigh angle).

### AVAILABLE COVERS

Vicair 02 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.

### PRODUCT PERFORMANCE

#### Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲▲
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

#### Positioning

Lateral stability	▲▲▲▲▲
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲▲

#### Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

#### Functionality & Application

Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲▲△
Available standard sizes	▲▲▲▲▲
Maximum user weight	▲▲▲▲△

### TECHNICAL SPECIFICATIONS

Thickness	6 and 10 cm
Average weight	725 grams
Max. user weight	250 kg
Machine washable	60°C max.

### CUSHION SIZES

Standard (Sizes >53 cm thickness 10cm only)

On request

		width											
Inch		14	15	16	17	18	19	20	21	22	23	24	
	CM	35	38	40	43	45	48	50	53	55	58	60	
depth	14	35											
	15	38											
	16	40											
	17	43											
	18	45											
	19	48											
	20	50											
	21	53											
	22	55											
	23	58											
	24	60											

070069-OVR-EN-REV01/2024