

 Breezy®



ה
ר
ר
ח
ה
ה
ה

ADJUSTMENTS MADE EASY



Easy to operate sliding seat for safe and easy transfer in and out of the wheelchair. By sliding the seat backwards, the user is positioned correctly in the wheelchair reducing the need for heavy lifting.



Individual adjustments and resetting of the wheelchair are simplified through the clear measuring scales.



Easy to use and identify handles and grips make adjusting the chair simple.

FLEXIBILITY AND QUALITY

Infinite adjustment possibilities built in....

From angle adjustable push bars, a sliding seat, easy to operate tilt and detachable swing-away legrests, the **Cirrus G5** is designed to provide an infinite amount of possibilities through a wide range of built in options and accessories.

Quality has been a key focus in the development of Cirrus G5. This includes high quality materials and technical solutions to provide a high quality wheelchair for users, carers and therapists.

Cirrus G5 meets the users need for comfort and stability, the therapists need for making individual adaptations and the technicians need for a high quality and flexible product.

COMFORT IS KEY

Excellent Comfort and Stability

The **Cirrus G5** provides high levels of comfort combined with an infinite level of adjustment possibilities including comfortable upholstery and arm rests that follow the back when reclining, the Cirrus G5 provides optimal support and comfort all day long.

Cirrus G5 has numerous adjustment options that ensure a good and stable seat position. The back recline and tilt provides an easy and comfortable way to change position throughout the day. These adjustments can be done whilst the user is sitting in the wheelchair and doesn't require the use of any tools.



CIRRUS G5

Designed to provide comfort and flexibility

throughout the day - both in activity and in rest.

ARMRESTS

Armrests follow the backrest angle to provide support & comfort in all positions. The armrest can be flipped up and removed to make transfer and transport easier.

BACKREST

The angle adjustable & individually adaptive back frame allows for a wide range of seating options to be fitted.

The standard comfort back has extra soft upholstery to ensure continued comfort. (pictured)

HEADREST OPTION

Comfort headrest designed with easy to adjust built-in straps ensures good support, both laterally and in-depth.

SEATING

A sliding seat pan allows for a wide range of seating cushions to be fitted easily.

A standard contoured seat cushion and pelvic seat plate provides pressure relief, stability and comfort for users. (pictured)

PUSH BAR

Angle adjustable and ergonomically designed. Standard on all Models.

LEGRESTS

The Cirrus G5 comes standard with angle adjustable legrests which follow the knees natural centre of rotation. Ankle adjustments can be done without tools.

TILTING

Double linkage tilting ensures safe and comfortable tilting down to 33°. This combined with the 35° back angle provides varied seating position options.





SPECIFICATIONS

Seat Width	42, 48, 54 cm
Seat Depth	39-55 cm
Seat Height (without cushion)	38 - 55 cm
Back Height	55 and 65 cm
Back Angle	0° - 35°
Tilt	0° - 33°
Armrest Height	22 - 33 cm
Total width	Seat width + 22 cm (+24 with drum brake)
Weight of chair	from 40 kg
Max user weight	135 kg