

**PERFECT TO USE IN COMBINATION WITH
LIBERTY BACK AND ANATOMIC BACK**

Hook and loop strips to
attach cushion

Adjustable straps
to adjust to back
contours. Number of
straps depends on
Strap back length

Easy and secure full length hook and
loop attachment to wheelchair canes,
with screw fixation for additional
security

Double zip at the front for addition
of lateral support

Double zip at the back to allow
fixation of stabilisation bar



- ▲ **Easy to apply on
manual wheelchairs**
- ▲ **Perfect positioning
and comfort qualities**
- ▲ **Very lightweight**

The Vicair Strap back is a perfect base for Vicair wheelchair back cushions. Quickly installed and combined with your choice of Vicair back cushion it is a high quality alternative to your standard wheelchair back system. The adjustable straps enhance the function of the Vicair back cushions, by offering thoracic stability and backrest positioning.

USED WITH

Vicair Liberty back, Vicair Anatomic back, Vicair Multifunctional O2 back (in combination with the Vicair Back Connector)

TECHNICAL SPECIFICATIONS

Max. user weight	150 kg
------------------	--------

rehab**hire** & sales

320 Lorimer Street Port Melbourne Victoria 3207

t: 03 9646 7200 | e: contact@rehabhire.com.au | rehabhire.com.au

070069-SB-EN-REV01/2023