



aspire

BATHROOM RANGE

aspire for...

Safety

Comfort

Support

Shower Chairs

Shower Stools

Over Toilet Aids

Toilet Surround Frame

Raised Toilet Seats



FEATURING
THE NEW
MAXI
RANGE

www.aspirecare.com.au
1300 133 120





- DVA Approved
- AS/NZ Approved
- Shower Chairs are registered in the ARTG, AUST R ARTG 230159/230165



Shower Chair - Wide

MODELS

BTS118000	Shower Chair with Arms - Steel
BTS118005	Shower Chair with Arms - Aluminium
BTS118010	Shower Chair with Arms - Wide - Steel
BTS118015	Shower Chair with Arms - Wide - Aluminium



1	Safe Working Load (S.W.L) of 200kg. Tested to AS/NZS3973.7:2009
2	Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
3	Height adjustable legs - Easy push button height adjustment for changing needs
4	Additional height increments on one leg - Accommodates uneven surfaces

5	All frames are individually treated for corrosion resistance
6	Raised armrests - Assist the user transferring on and off the chair
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
8	Heavy duty, reinforced frame accommodates very high SWL

SIZE	OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	UNIT WEIGHT	S.W.L
470mm	735-910 mm	512 mm	570 mm	400 mm	440 mm	470 mm	375-550 mm	Alum: 3.1 kg Steel: 3.4 kg	200 kg
525mm (wide)	735-910 mm	570 mm	570 mm	400 mm	440 mm	525 mm	375-550 mm	Alum: 3.3 kg Steel: 3.9 kg	200 kg



- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ MAXI is registered in the ARTG, AUST R ARTG 306118



2



4



7



3



5

MODEL


BTS118020 | MAXI Shower Chair

5 Year Warranty
Corrosion

1 Year Warranty
Structure

1	Safe Working Load (S.W.L) of 310kg.
2	Super heavy duty width adjustable frame to 775mm
3	Soft padding on seat and arms for postural support and comfort
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant

5	Height adjustable legs with easy circlip adjustment points for changing needs
6	Raised armrests – Assist the user transferring on and off the chair
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability





- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ Shower Stools are registered in the ARTG, AUST R ARTG 230166/230167



BTS119025





MODELS

- BTS119000 Shower Stool - Steel
- BTS119005 Shower Stool - Aluminium
- BTS119020 Shower Stool with Padded Seat (PS) - Steel
- BTS119025 Shower Stool with Padded Seat & Back (PS & B) - Steel



1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Clip on seat or padded back attachment – For easy assembly and cleaning
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	All frames are individually treated for corrosion resistance

5	Seat with anterior tilt – Facilitates standing from sitting and provides greater comfort for stiff or painful hips
6	Raised armrests – Assist the user transferring on and off the stool
7	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces

OVERALL HEIGHT	SEAT HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	UNIT WEIGHT	S.W.L
								
670-855 mm	470-655 mm	545 mm	400 mm	310 mm	375 mm	495 mm	Alum: 3.1 kg Steel: 3.8 kg	175 kg



- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ Shower Chairs are registered in the ARTG, AUST R ARTG 230166/230167



MODELS

- BTS119010 | Shower Stool - Wide - Steel
- BTS119015 | Shower Stool - Wide - Aluminium

5
Year
Warranty
Corrosion

1
Year
Warranty
Structure

1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Clip on seat attachment – For easy assembly and cleaning
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	All frames are individually treated for corrosion resistance
5	Seat with anterior tilt – Facilitates standing from sitting and provides greater comfort for stiff or painful hips

6	Raised armrests - Assist the user transferring on and off the stool
7	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces
9	Comfortable wide seat with drainage holes – For increased support and hygiene
10	Reinforced wide frame – Provides stability for larger users

OVERALL HEIGHT	SEAT HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	UNIT WEIGHT	S.W.L
 575-765 mm	 425-615 mm	 570 mm	 500 mm	 365 mm	 420 mm	 525 mm	Alum: 3.1 kg Steel: 3.8 kg	 175 kg



- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ MAXI is registered in the ARTG, AUST R ARTG 306118



2



4



6



5



7

MODEL


BTS119030 | MAXI Shower Stool

5 Year Warranty
Corrosion

1 Year Warranty
Structure

1	Safe Working Load (S.W.L) of 310kg.
2	Super heavy duty width adjustable frame to 775mm
3	Soft padding on seat and arms for postural support and comfort
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant

5	Height adjustable legs with easy circlip adjustment points for changing needs
6	Raised armrests – Assist the user transferring on and off the stool
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability





- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ Toilet Surround Frame is registered in the ARTG, AUSTRALIA ARTG 230170



MODEL

BTT155100 | Toilet Surround - Aluminium

1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Height adjustable frame – Easy push button height adjustment for changing needs
3	Raised armrests – Assist the user transferring on and off the toilet

4	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
5	Additional height increments on one leg – Accommodates uneven surfaces

OVERALL HEIGHT	OVERALL DEPTH	OVERALL WIDTH	WIDTH BETWEEN ARMS	UNIT WEIGHT	S.W.L
					
535 - 660mm	435mm	500mm	450mm	3 kg	175 kg



- DVA Approved
- AS/NZ Approved
- Over Toilet Aids are registered in the ARTG, AUST R ARTG 230169/230170

Over toilet Aid - Wide

MODELS

- BTT145000 Over Toilet Aid - Steel
- BTT145005 Over Toilet Aid - Aluminium
- BTT145010 Over Toilet Aid - Wide - Steel
- BTT145015 Over Toilet Aid - Wide - Aluminium

5
Year
Warranty
Corrosion

1
Year
Warranty
Structure

1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Raised armrests – Assist the user transferring on and off the over toilet aid
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	Comfortable, contoured seat with clip on attachment – For easy assembly and cleaning. Padded seat available, sold separately (BTT145110)

5	All frames are individually treated for corrosion resistance
6	Removable, moulded splash guard provided – Easily removed for transport and cleaning
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces

SIZE	OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	UNIT WEIGHT	S.W.L
460 mm	700-880 mm	512 mm	625 mm	435 mm	370 mm	460 mm	450-630 mm	Alum: 3.4 kg Steel: 4.2 kg	175 kg
525 mm (Wide)	700-880 mm	570 mm	625 mm	435 mm	370 mm	525 mm	450-630 mm	Alum: 3.4 kg Steel: 4.2 kg	175 kg



- ✓ DVA Approved
- ✓ AS/NZ Approved
- ✓ MAXI is registered in the ARTG, AUST R ARTG 306118


5
Year
Warranty
Corrosion

1
Year
Warranty
Structure

MODEL
BTT145020 | MAXI Over Toilet Aid

1	Safe Working Load (S.W.L) of 310kg.
2	Super heavy duty width adjustable frame to 775mm
3	Flip up padded seat for ease of cleaning
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant

5	Height adjustable legs with easy circlip adjustment points for changing needs
6	Raised armrests – Assist the user transferring on and off the over toilet aid
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability





✓ DVA Approved

✓ AS/NZ Approved

✓ Toilet Seat Raisers are registered in the ARTG, AUST R ARTG 242921



MODELS

BTT149800	Toilet Seat Raiser - 50mm - No Lid
BTT150200	Toilet Seat Raiser - 100mm - No Lid
BTT150300	Toilet Seat Raiser - 150mm - No Lid
BTT149900	Toilet Seat Raiser - 50mm - With Lid
BTT150100	Toilet Seat Raiser - 100mm - With Lid
BTT150400	Toilet Seat Raiser - 150mm - With Lid




1	Three heights available – helps facilitate sit-to-stand transfers for safe and comfortable toileting
2	Contoured plastic seat – Comfortable and durable
3	Available with or without lid

4	Fully sealed structure – Easy to clean for increased infection control
5	Generous cut outs at the front and rear – Facilitates personal care and hygiene
6	Two adjustable mounting brackets – Provide an adjustable, secure fit

OVERALL HEIGHT	UNIT WEIGHT	S.W.L
 2" (50mm) 4" (100mm) 6" (150 mm)	 1.4 kg	 200 kg



SHOWER COMMODES



BATH SPONGES



BATH LIFTS



BATH/SHOWER TRANSFER BENCHES



BATH STEPS



GRAB RAILS



HAND HELD SHOWERS



SWIVEL BATHERS



BATH BOARDS

CRUTCHES, FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



RAMPS



AIR MATTRESSES & FOAM MATTRESSES



COMMODES, URINALS & BED PANS



LIFTERS & SLINGS



BEDS & BED POLES



DAILY LIVING AIDS



SHOWER COMMODES



BATH BOARDS & BENCHES



SEAT WALKERS



POWER WHEELCHAIRS




PRESSURE CUSHIONS



LIFT & RECLINE CHAIRS



Distributed by:



ASBR26.26 300719

