

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

PRIDE[®]
MOBILITY

PATHRIDER[®] 130 XL

Mobility Scooter



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Unstoppable Performance

The Pathrider® 130 XL is designed for those seeking a reliable ride with the power to handle uneven terrain. Equipped with a robust 4-pole motor, this scooter effortlessly navigates both paved paths and rugged surfaces, reaching speeds of up to 10km/h. The high-back reclining seat and wraparound delta tiller ensure smooth and simple steering, while front and rear suspension provide a stable, comfortable ride over outdoor terrain.

Specifications

Width ³	636mm
Length ³	1360mm
Seat Width	495mm
Seat Depth	445mm
Seat to Deck Height	470mm suspension seat post
Seat to Ground Height	567mm suspension seat post
Ground Clearance ⁵	90mm at transaxle
Turning Radius	1700mm
Tyres Front	330.2mm
Tyres Rear	330.2mm
Range per Charge ¹	<25km
Maximum Speed ¹	10kph
Battery Charger	8A off-board
Battery Requirements ⁴	55AH x 2
Battery Weight	16.8kg Each
Weight without Battery ²	99.8kg
Colours	Candy Apple Red Mercury Silver
SWL	158kg
Compliance	AS3695.2:2013
Warranty	2 years Batteries: 12 months



Full Suspension



Easy Drive Tiller



Pneumatic Tyres

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Width of base (not including seat).
4. Deep cycle (AGM or gel cell type recommended)
5. Does not account for anti-tip clearance.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

rehab**hire** & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire** & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider