

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

**PRIDE**<sup>®</sup>  
MOBILITY

# **PATHRIDER<sup>®</sup> 140 XL**

Mobility Scooter



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## Built for Adventure!

The Pathrider® 140 XL combines power and precision in an exhilarating ride. Equipped with a robust 24-volt DC motor, full suspension, hydraulic-sealed brake system, and large pneumatic tyres, this scooter is ready to conquer any terrain. Premium features like the adjustable wraparound tiller, high-back pillow-top seat, and bright LED lighting add a touch of luxury. Whether cruising down the street or exploring the outdoors with family and friends, the Pathrider® 140 XL is built for the ultimate outdoor adventure.

## Specifications

Width <sup>3</sup>	692mm
Length <sup>3</sup>	1428mm
Seat Width	495mm
Seat Depth	445mm
Seat to Deck Height	470mm
Seat to Ground Height	667mm – 705mm
Ground Clearance <sup>5</sup>	133mm at midframe
Turning Radius	Right: 2100mm   Left: 2080mm
Tyres   Front	330.2mm
Tyres   Rear	330.2mm
Range per Charge <sup>1</sup>	<35km
Maximum Speed <sup>1</sup>	10kph
Battery Charger	8A off-board
Battery Requirements <sup>4</sup>	70AH x 2
Battery Weight	23.5kg Each
Weight without Battery <sup>2</sup>	117.7kg
Colours	Candy Apple Red Silver   Yellow
SWL	181kg
Compliance	AS3695.2:2013
Warranty	2 years   Batteries: 12 months



133mm Ground Clearance



Full Suspension



High Back Pillow Top Seat

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Width of base (not including seat).
4. Deep cycle (AGM or gel cell type recommended)
5. Does not account for anti-tip clearance.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.



# rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider