


ottobock.

The new Start product range presents itself

Lighter, more functional, more versatile







The new generation of our proven Start series presents itself

The demands users make on their technical devices are as complex as their individual needs and everyday habits. Taking this into account is the goal of the new generation of our proven Start wheelchairs with numerous updates. Expert feedback and patient experience inspired us to make targeted adaptations that offer a wide range of users exactly the support they need. You can take a detailed look at our special optimisations on the following pages.



Antibacterial effect of the push handles

All models of the Start series are equipped with antibacterial push handles. The antibacterial effect of the push handles has been verified by ACTEGA DS GmbH. In this process it was possible to demonstrate under laboratory conditions that the push handles, which are produced from the material ProvaMed®, are germ-resistant and thus have an antibacterial effect. The antibacterial effect represents an important property especially in institutional domains for preventing transfer of bacteria.

Start M2 - Start 4 configurable model

The Start M2 is a fully configurable, high-quality wheelchair model for which all of the options of the Start range are available. This model offers the maximum number of combination possibilities with regard to dimensions, colours and adjustments which makes it possible to provide even large and very large users with an optimal device. Drum brakes are also available options, as are elevating footrests or amputation legrests. The Start M2 also has an adaptive component. The rear wheel adapter allows you to adjust the wheelchair precisely to the user's needs.

There are many options available with the new Start M2 for the fitting of stroke patients with hemiplegia and cerebral palsy patients. The Start M2 is just the right choice for custom fittings in the lightweight segment.



Technical data

User Weight	130 kg
Seat width	380 / 405 / 430 / 455 / 480 / 505 mm
Seat depth frame size 1	405 / 430 / 455 mm
Seat depth frame size 2	455 / 480 / 505 mm
Back height	400 / 425 / 450 / 475 mm
Anterior seat height	380 – 510 mm
Posterior seat height	380 – 510 mm
Seat angle	selectable
Lower leg length	250 – 470 mm (continuously adjustable)

Product highlights

- Fully configurable lightweight wheelchair
- Patented, tool-free, quick seat depth adjustment (up to 50 mm)
- Very low weight with full Ottobock quality (from 14.7 kg total weight, from 9.6 kg transport weight*)
- Antibacterial push handles**
- Excellent component compatibility
- New adjustable economy back upholstery
- Larger centre of gravity adjustment range
- Larger selection of seat depths, including for people with a height over 2 m
- Larger selection of seat widths
- More options for height adjustment and seat tilt
- Wider foot space for hemiplegics, which allows “foot propelling”
- New one-hand drive for hemiplegics
- Numerous updated and improved accessories
- Two different frame lengths for optimal seat depth adjustment
- Selection of 6 different colours (anthracite metallic, jet black RAL 9005, ice blue, silver metallic, saturn blue metallic, white RAL 9016)
- Side panel with armrest, height can be adjusted (without tools) or set (with tools), adjustable depth

* The transport weight is the total weight of the wheelchair minus the weight of the rear wheels, footrests and seat cushion.

** The antibacterial effect of the push handles has been verified by an independent institute.

Start M2S - Start 4 standard and Start 4 adjustable model

In order to ensure the rapid availability of our top model Start M2, there are many stock versions, which we refer to as M2S. The “S” stands in this case for “stock”. For you this means that these models are available quickly thanks to the optimum logistics chain and can thus be used quickly by patients.

Technical data



	M2S Standard	M2S Adjustable
User weight	130 kg	130 kg
Seat width :	405 / 455 / 505 mm	380 / 405 / 430 / 455 / 480mm
Seat depth :	405 / 430 / 455 mm	405 / 430 / 455 mm
Back height :	400 / 425 / 450 / 475 mm	400 / 425 / 450 / 475 mm
Anterior seat height	430 - 520 mm	430 - 520 mm
Posterior seat height	410 - 505 mm	410 - 505 mm
Seat angle	0 – 5°	0 – 5°
Lower leg length	280 – 500 mm (continuously adjustable)	280 – 500 mm (continuously adjustable)
Rear wheels	24" standard	24" standard
Caster wheels	8" PU (200x50mm)	6" PU
Side Panel with armrest	Height adjustable (with tools)	Height adjustable (10cm)(push button)
Frame colour	Silver Metallic	Silver Metallic
Back	With Lumbar curve	Straight
Back upholstery	Standard	Adjustable economy
Push handles	Standard	Standard
Seat upholstery	Standard	Standard

Product highlights

- Stock version of the fully configurable Start M2
- Patented, tool-free, quick seat depth adjustment (50 mm)
- Very low weight with full Ottobock quality (from 14.7 kg total weight, from 9.6 kg transport weight*)
- Antibacterial push handles**
- Excellent component compatibility
- New adjustable economy back upholstery
- Short delivery times for stock versions (quick availability)
- Aluminium frame
- Frame colour: silver metallic
- Side panel with armrest, height can be adjusted (without tools) or set (with tools) and adjustable depths

* The transport weight is the total weight of the wheelchair minus the weight of the rear wheels, footrests and seat cushion.

** The antibacterial effect of the push handles has been verified by an independent institute.

Start M2 options



Side panel with height-adjustable armrest



One-Arm drive



Armrest with swivel unit



Elevating footrest



Amputation Legrest



Add-on drives (i.e the new eSupport)

Otto Bock Australia

Suite 1.01, 62 Norwest Blvd, Baulkham Hills NSW 2153

www.ottobock.com.au - healthcare@ottobock.com.au

rehab**hire** & **sales**

320 Lorimer Street Port Melbourne Victoria 3207

t: 03 9646 7200 | e: contact@rehabhire.com.au | www.rehabhire.com.au