

## ULTRA-LIGHTWEIGHT, EASY TO USE AND COMFORTABLE

**Compartment barrier** ..... gently limits sliding forward

### Equal pressure distribution on the seating surface



## Auto-adjusting valve

- ▲ **Air valve automatically adapts the cushion to your body shape and weight**
- ▲ **Two compartment layout limits sliding**
- ▲ **Easy to clean and disinfect**

The Vicar Liberty is an easy to use and comfortable wheelchair cushion. It is designed for users who require comfort in combination with a low to medium level of support and pressure injury prevention. The cushion is ready-for-use, rain proof and with two compartments it provides an even pressure distribution. This makes it an ideal cushion in a large variety of situations.

The automatic valve makes sure that the cushion adjust to every seating plane and every body contour.

This Vicair cushion is extremely lightweight.

## PRODUCT PERFORMANCE

### Skin protection / Anti decubitus

Pressure redistribution	▲▲▲△△
Shear force reduction	▲▲△△△
Microclimate control	▲▲△△△
Impact & vibration damping	▲▲▲△△

## Positioning

Lateral stability	▲▲▲△△
Frontal stability	▲▲▲▲△
Sliding reduction	▲▲▲▲△

Allows manual adjustments	N/A
---------------------------	-----

Increasing endurance ▲▲▲△△

## Maintenance & Safety

Durability	▲▲▲▲△
Reliability	▲▲▲▲△
Functionality when misused	▲▲▲▲△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

### Functionality & Application

Comfort	▲▲▲▲△
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲△△
Maximum user weight	▲▲▲△△

## TECHNICAL SPECIFICATIONS

Thickness	5 cm
Average weight	350 grams
Max. user weight	150 kg

## CUSHION SIZES

- Standard

■ Not available

[illegible]