

JAY FUSION™

The JAY® Fusion™ is an adjustable cushion designed for skin protection that features JAY® Flow™ fluid or JAY® air insert options in a package that's comfortable, stable and supportive.

Common Client Characteristics

- May be at high risk of skin breakdown
- May not be able to perform independent weight shifts
- May have changing postural needs over time
- May have aggressive positioning needs
- May be in a wheelchair for extended periods of time
- May experience weight loss



Fusion Rating Scale

Side-to-Side Stability	4.5
Forward Stability	4
Pressure Management	4.5

Product Code				
CUSHION with fluid	CUSHION with single air	CUSHION with dual air	COVER	BASE
JFUSIONXXXX	JFUSIONXXXXAS	JFUSIONXXXXAD	CJFUSIONXXXX / CJFUSIONXXXXA	BJFUSIONXXXX
INNER COVER	AIR BLADDER SINGLE VALVE	AIR BLADDER DUAL VALVE	FLUID	
ICJFUSIONXX	ASAJFUSION / ASBJFUSION / ASCJFUSION	ADSJFUSION / ADBJFUSION / ADCJFUSION	FAJFUSION / FBJFUSION / FCJFUSION	
Component	Feature	Design Intent	Clinical Value	
BASE	Pre-contoured closed cell foam base	Provides lateral and forward stability that is easy to modify	Encourages orthopedic alignment, increases sitting tolerance and accommodates changing user needs	
	Pre-contoured, soft foam layer with posterior pelvic wall	Provides positioning of the pelvis and thighs and reduces fluid migration	Postural support with comfort	
INSERT	Optional JAY® Flow™ fluid pad	Consistently ensures proper fluid placement by reducing fluid migration	Reduces risk of skin breakdown	
	Optional JAY® AIR INSERT air bladder available in single valve or dual valve options	Allows personalisation of air levels to optimize immersion and accommodate asymmetric postures	Reduces risk of skin breakdown	
INNER COVER	Inner cover with Aquaquad™ zipper and anti-wicking thread	Protects foam base from moisture build up and is easy to clean	Helps with incontinence management	
COVER	Ionic silver thread stretch outer cover	Reduces surface tension and resist bacterial growth	Helps prevent bacterial growth and reduces risk of skin breakdown	

Technical Specifications

Waterlow	25+
Width Range	14" to 24"
Depth Range	14" to 24"
Weight	2.4kgs 16"x18" w/fluid)
Height	4"
	3" Reduced Profile
Weight Capacity	150kgs (14-21" widths)
	227kgs (22-24" widths)
JAY Your Way Modifications Avail.	Yes
Positioning	Yes

	WIDTH													
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
DEPTH	12"	N/A	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A	N/A	N/A	
	13"	N/A	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	
	14"	N/A	JYW	●	●	●	●	●	JYW	JYW	JYW	JYW	N/A	N/A
	15"	N/A	JYW	●	●	●	●	●	●	●	JYW	JYW	JYW	JYW
	16"	N/A	JYW	●	●	●	●	●	●	●	JYW	JYW	JYW	JYW
	17"	N/A	JYW	●	●	●	●	●	●	●	●	●	JYW	JYW
	18"	N/A	JYW	●	●	●	●	●	●	●	●	●	JYW	JYW
	19"	N/A	JYW	JYW	●	●	●	●	●	●	●	●	●	●
	20"	N/A	JYW	JYW	●	●	●	●	●	●	●	●	●	●
	21"	N/A	JYW	JYW	●	●	●	●	●	●	●	●	●	●
	22"	N/A	N/A	N/A	JYW	●	●	●	●	●	●	●	●	●
	23"	N/A	N/A	N/A	N/A	JYW	●	●	●	●	●	●	●	●
24 "	N/A	N/A	N/A	N/A	N/A	JYW	●	●	●	●	●	●	●	

● Standard size available JYW Size available through JAY Your Way N/A Size not available