

Feel in control with Comet Alpine+

For those who enjoy adventures across terrains too challenging for a traditional scooter, the Comet® Alpine+ offers the perfect solution.

Its robust and powerful driving characteristics gives you the freedom to discover the great outdoors with minimal effort in continuous comfort and safety.

Up to
10 km/h
max speed

Up to
58 km
max distance
on a single
charge

160 kg
max user
weight



Premium performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension helps enable a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance*
- Anti-splash guards



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to help prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Adjustable steering column
for responsive control
and maximum comfort



Seat suspension module



Seat riser



LCD dashboard for clear, easy to read information.



Add to your scooter

Tailor the Invacare Comet Alpine to your taste and needs. Turn over to explore the optional extras available.