

Invacare Scooters

Feel the Freedom



Yes, you can.®

Content

Free 24 Hour Roadside Assistance	3
Colibri	4
Leo	6
Pegasus Metro	8
Pegasus Pro	10
Comet Pro (NZ Only)	12
Comet Alpine+	14
Comet Ultra	16
Accessories	18



Keeping you moving
Free 24 Hour
Roadside Assist



New Zealand



AA help is at hand 365 days of the year

Invacare has joined up with the people at AA Roadservice cover to bring you Peace of mind. With any Invacare mobility scooter purchase, you will receive 24 hour, 365 days per year, roadside assistance, for three years, FREE! So not only can you enjoy all the great features offered by our range of mobility scooters, you can rest assured knowing that if anything goes wrong with your mobility scooter, help is just an 0800 number away.

How do I become eligible for the 3 years FREE offer?

Simple. If you purchase any Invacare mobility scooter, you need to ensure your Invacare approved dealer registers for your warranty online at the time of purchase*. Once your dealer has registered for your warranty, you will be automatically signed up for the program and can call 0800 500 222 any time you need to!



Australia



24/7 Road Service Australia

Invacare offers 1 year free breakdown assistance when you purchase your new scooter. To qualify for roadside assistance customers need to have current and active policies with 24/7 Road Service Australia.

Roadside assistance provides assistance only when:

- Customer has active policy
- Broken down (Mechanical Breakdown ONLY) outside of home

Roadside assistance provides a towing service to bring the scooter home for the customer. We do not offer repairs, we do not do maintenance nor offer assistance of any kind when customer are located in their home.

To be eligible you need to ensure your Invacare approved dealer registers your warranty online at the time of purchase. Once your dealer has registered your warranty, you will be automatically signed up for cover and can call 1300 550 279 any time you need to!

For further information ring 24/7 Road Services Australia on 1300 550 279.



Feel in control with Invacare Colibri®

Introducing the Invacare Colibri®, our lightest and most transportable scooter. Designed to complement your style with a choice of colours and improve your independence.

Compact and easy to drive, this scooter has been designed with a patented Litelock™ system making it easy to transport into the back of your car. The swivel seat makes getting on and off your scooter simple. The Colibri® is

incredibly simple to start and steer. It has been designed to give you a comfortable and safe drive with the ability to manoeuvre in tight spaces.



Easy transportation with the Invacare Litelock™ system



Easy to use and transport

- Travel easily in style
- Ergonomic steering for added control
- Simple to steer and drive




Safety in mind

- Reflectors
- Additional optional extras available



Swivel seat for easy transfer



Soft, padded seat for enhanced comfort




Adjustable steering for improved seating position



Our **lightest** and
most transportable
scooter



Adjustable steering
for greater comfort



Easy battery charge
in situ or remove
from scooter



Add to your scooter

Tailor the Invacare Colibri® to your taste and needs. Turn over to explore the optional extras available.

Feel in control with Invacare Leo®

The Invacare Leo® is a compact scooter that gives a smooth, safe drive. Simply dismantle it for easy transportation and take it with you for lots of exciting days out and those moments that matter.

Leo® is packed with many safety and performance features. Its intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. A visual warning for the battery level display, alerts

with an audible warning (beep) when the battery needs recharging. Incredibly easy to use the Leo® puts you in control with confidence.



Easy to use and transport

- Dismantles into five parts for easy transportation
- Ergonomic steering for added control
- Simple to steer and drive



Safety in mind

- Two-step disengaging lever to prevent free-wheeling
- Brake lights, headlights and indicators
- Additional optional extras available



Anti-splash guards

protects the scooter from dirt and water



Adjustable steering for greater comfort



Easy battery charge

with off-board charger kit



Swivel seat for easy transfer



Easy to transport in your car



Two-step disengaging lever to prevent free-wheeling



Add to your scooter

Tailor the Invacare Leo® to your taste and needs. Turn over to explore the optional extras available.

Feel in control with Pegasus Metro

Introducing our Pegasus Metro, a high-quality scooter designed to enhance safety, comfort and added driving performance. The Pegasus Metro offers everything you need from a mobility scooter.

Compact and easy to drive, this scooter handles a wide variety of surfaces. Designed with a powerful motor, the Pegasus Metro drives across uneven ground, around corners and over most obstacles without compromising

safety or comfort. Incredibly simple to steer and start, the Pegasus Metro puts you in total control.



Advanced driving control

- Ergonomic steering for added control
- Designed for all-terrain performance*
- Anti-splash guards



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to help prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Swivel seat for easy transfer



High quality lighting



Two-step disengaging lever to help prevent free-wheeling



Designed to help **provide traction around corners** and safe, all-terrain performance*



Adjustable steering for maximum comfort



Add to your scooter

Tailor the Invacare Pegasus Metro to your taste and needs. Turn over to explore the optional extras available.

* Please refer to the User Manual for limitations of use.

Feel in control with Pegasus Pro

Introducing our Pegasus Pro, a high performing scooter designed for your comfort, safety and reliability wherever you go.

Ideal for longer journeys, the Pegasus Pro features advanced suspension and large wheels to ensure you drive in comfort and safety even across uneven ground,

around corners or over obstacles. Incredibly simple to steer and start, the Pegasus Pro puts you in total control.

up to
15 km/h
max speed

up to
52 km
max distance
on a single
charge

160 kg
max user
weight



Advanced driving control

- Ergonomic steering for added control
- Designed for all-terrain performance*
- Anti-splash guards



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever helps prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Swivel seat for easy transfer



High quality lighting



Two-step disengaging lever to help prevent free-wheeling



Designed to help **provide traction around corners** and safe, all-terrain performance*



Adjustable steering for maximum comfort



Add to your scooter

Tailor the Invacare Pegasus Pro to your taste and needs. Turn over to explore the optional extras available.

Feel in control with Comet Pro

ONLY AVAILABLE
IN NEW ZEALAND

The Comet Pro is our most advanced outdoor scooter. Designed with exceptional features and options to take you wherever you want to go in comfort, safety and style. It's a design to provide drive performance and help enable control without compromise.

Comet Pro features exceptional driving features and modern aesthetics to deliver premium performance and control. The 2 x 12 V / 75 Ah batteries allow users to enjoy challenging journeys beyond the city's limits.

It also has large front and rear wheels, combined with an advanced suspension system that helps enable safe negotiation of curbs and obstacles smoothly and comfortably.

Up to
15 km/h
max speed

Up to
55 km
max distance
on a single
charge

160 kg
max user
weight



Premium performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension helps enable a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance*
- Anti-splash guards



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to help prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Exceptional suspension



Reinforced seat column



Height and width adjustable armrest



Auto-speed reduction

helps reduce speed around corners for added safety



Adjustable steering column
for responsive control and maximum comfort



Add to your scooter

Tailor the Invacare Comet Pro to your taste and needs. Turn over to explore the optional extras available.

* Please refer to the User Manual for limitations of use.

Feel in control with Comet Alpine+

For those who enjoy adventures across terrains too challenging for a traditional scooter, the Comet® Alpine+ offers the perfect solution.

Its robust and powerful driving characteristics gives you the freedom to discover the great outdoors with minimal effort in continuous comfort and safety.



Premium performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension helps enable a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance*
- Anti-splash guards



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to help prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Adjustable steering column
for responsive control and maximum comfort



Seat suspension module



Seat riser



LCD dashboard for clear, easy to read information.



Add to your scooter

Tailor the Invacare Comet Alpine to your taste and needs. Turn over to explore the optional extras available.

Feel in control with Comet Ultra

Ideal for those with larger body shapes, the Comet® Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality.

Packed with all of the safety features and driving characteristics seen throughout the Comet® family, the Comet® Ultra will get you to

your destination quickly and efficiently, in complete comfort and control.

Up to
10 km/h
max speed

Up to
49 km
max distance
on a single
charge

220 kg
max user
weight



Improved performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension ensures a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



Auto-speed reduction

helps reduce speed around corners for added safety



Safety in mind

- Automatic speed reduction for safe cornering
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Improved suspension



Reinforced seat column



Height and width adjustable reinforced armrest for extra support



Adjustable steering column

for responsive control and improved comfort



Add to your scooter

Tailor the Invacare Comet Ultra to your taste and needs. Turn over to explore the optional extras available.

Driving Comfort

Part No.	Description	Comet	Pegasus	Leo	Colibri
Seat suspension					
A seat suspension kit adds extra comfort during driving.					
Suspension Height					
SP1577212	Low	425mm up to 100kg			●
SP1577213	Medium	450mm up to 100kg			●
SP1577214	High	475mm up to 100kg			●
1522456	Low		●		
1522457	Medium		●		
1522458	High		●		
SP1541627	Low			●	
SP1541628	Medium			●	
SP1541629	High			●	
SP1528986	Low	●			
SP1528987	Medium	●			
SP1528988	High	●			
1522459	Low	●			
1522460	Medium	●			
1522461	High	●			
375000-68600	Sun canopy Add the sun shade when needed.	●	●	●	
SP1543160	Wind screen Offers protection from the wind.	●	●	●	
1505731	Storage Bag This spacious bag offers extra storage room.		●		

Storage

Part No.	Description	Comet	Pegasus	Leo	Colibri
1637910	Lockable storage box Keeps your valuables secure while driving; box can be mounted in the front or rear of the scooter. (SP1552537 bracket required to fit storage box to the rear. No bracket required for front mounting)	●	●	●	
1637913	Rear lockable storage box (mounting bracket included)	●	●		
1505727	Rear basket Used to store your extra shopping.	●	●	●	●
SP1577430	Front basket (plastic)	●	●	●	●
1637915	XL Front basket (plastic)	●	●		

Maintenance

1509434	Storage cover			●	●
SP1518240	Protects your scooter from dirt and wet conditions.	●	●		

Carriers

SP1553851	Rollator carrier The rollator carrier is adjustable to carry your rollator safely.	●	●	●	
1505728	Walking frame carrier Keeps your walking frame stable during the drive.	●	●	●	
SP1521594	Cup holder Stores bottles or soft drinks; can be folded up when not in use.	●	●	●	
1572971	Oxygen bottle carrier			●	●
1505729	Securely stores your oxygen during driving.	●	●		
1528274	Cane and crutch holder Keeps your cane/crutch stable during driving. <small>*Must order SP1577455 when using with Colibri</small>	●	●	●	●
1677347	Mobile phone holder Keeps your phone hands-free while driving.	●	●	●	●

Accessories






Part No.	Description	Comet	Pegasus	Leo	Colibri
SP1541502	Extension Tube				
	Use the extension bracket when mounting accessories in the seat receiver.	●	●		
SP1577455	Accessory Bracket				
	Use the bracket when mounting accessories in the seat receiver.				●
SP1552537	Rear Bracket for Storage Box				
	Use the front basket for your shopping and mount the storage box behind your seat using this bracket.	●	●	●	
SB02	Single Bracket Accessory				
	2:1 Adapter that allows single stem accessories be fitted to a twin stem seat mount.	●	●	●	
1637917	USB Charger	●	●		
1677347	Mobile Phone Holder	●	●	●	●

Safety Options



1509362	Seat belt			●	●
SP1518274	Keeps you well positioned in the seat at all times.	●	●		
	Steel bumpers				
	Protect the Scooter shrouds from scratches and bumps.				
1638965	Front bumper, Pegasus		●		
1638966	Rear bumper, Pegasus		●		
1528298	Front bumper, Comet	●			
1638967	Rear bumper, Comet	●			
SP1521592	Safety vest	●	●	●	●
	Enhances your visibility to others.				
SP1543153	Personal key				●
SP1543154	For extra safety choose a personal key cylinder for your scooter, available with 2 personal keys.	●	●	●	
SP1518241	Large mirrors	●	●	●	
	Improve your safety - can be mounted left or right.				
1637911	Deluxe mirrors	●	●		
	Larger mirrors for your Pegasus or Comet				

Colour Options:












Colibri:

		
Jade Green Code: 1577170	Sapphire Blue Code: 1577171	Jasper Red Code: 1577169
		
Marble White Code: 1577168	Sunstone Code: 1577172	

Leo: Complete Scooter

	
Onyx Blue Code: HS-528-BU4	Ruby Red Code: HS-528-RD2






Pegasus Metro, Pegasus Pro, Comet Pro, Comet Alpine +, Comet Ultra:






			
Jade Green Code: 1608822	Allum Blue Code: 1608823	Electric Blue Code: 1637418	Jasper Red Code: 1608832
			
Carmine Red Code: 1637417	Crystal Copper Code: 1637416	Luxor Beige Code: 1637415	Iconic Chrome Code: 1637419
			
Diamond Silver Code: 1608831	Marble White Code: 1608821	Canary Yellow Code: 1608833	






Whatever your final choice, enjoy the freedom and confidence of being able to get around comfortably and securely. Find out more about our Scooter Range at www.invacare.com.au or www.invacare.co.nz. Our Scooter Range offers you the perfect solution for every excursion or adventure. Imagine being able to take long journeys in complete safety and comfort.

Technical data

For more information about this product, including the user manual, please visit our website

						
	Backrest width ▼	Seat depth ▼	Seat to floor height ▼	Backrest height ▼	Driving unit width ▼	Total length ▼
Colibri	465 mm	400 mm	370 - 420 mm	385 mm	505 mm	1010 mm
Leo	470 mm	410 mm	460 - 560 mm	475 mm	590 mm	1220 mm
Pegasus Metro	510 mm	470 mm	Min. 440 mm Max. 510 mm	500 mm	660 mm	1270 mm
Pegasus Pro	510 mm	460 mm	Min. 440 mm Max. 510 mm	400 mm	650 mm	1320 mm
Comet Pro ¹	510 mm	470 mm	Min. 440 mm Max. 510 mm	500 mm	665 - 685 mm	1450 mm
Comet Alpine+	510 mm	460 mm	Min. 440 mm Max. 510 mm	630 mm	665 - 685 mm	1450 mm
Comet Ultra	660 mm	535 mm	Min. 440 mm Max. 510 mm	515 mm	665 - 685 mm	1450 mm

						
	Empty weight ▼	The heaviest part ▼	Max. User weight ▼	Battery capacity ▼	Turning diameter ▼	Max. Climbable obstacle height ▼
Colibri	49 kg	15 kg	136 kg	2 x 18 Ah	2200 mm	45 mm
Leo	83 kg	46 kg	136 kg	2 x 36 Ah	2620 mm	60 mm
Pegasus Metro	110 kg	57 kg	136 kg	2 x 12 V/40 Ah (C20) AGM	2600 mm	80 mm
Pegasus Pro	136 kg	62 kg	160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm
Comet Pro ¹	136 kg	63 kg	160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm
Comet Alpine+	143 kg	63 kg	160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm
Comet Ultra	148 kg	63 kg	220 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm

						
	Ground clearance ▼	Max. Safe slope ▼	Motor capacity ▼	Max. Driving distance ▼	Max. Driving speed ▼	Warranty ▼
Colibri	50 mm	6°	200 W	Up to 16 km	8 km/h	5 Years
Leo	70 mm	10°	240 W	38 km	8 km/h	5 Years
Pegasus Metro	80 mm	10°	550 W / 1500 W ²	54 km	10 km/h	5 Years
Pegasus Pro	100 mm	10°	550 W / 1600 W ²	52 km	10 km/h	5 Years
Comet Pro ¹	100 mm	10°	550 W / 1600 W ²	55 km	15 km/h	5 Years
Comet Alpine+	100 mm	10°	650 W / 1800 W ²	58 km	10 km/h	5 Years
Comet Ultra	100 mm	10°	650 W / 1800 W ²	49 km	10 km/h	5 Years

¹ NZ Only

² Peak performance

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au