

## THE OPTIMAL CUSHION FOR WHEELCHAIR USERS WITH BODY ASYMMETRY

**Follows the body's contours**  
for equal pressure distribution | support IT-coccyx area | managing pelvic asymmetry

**Pre-ischial bar limits sliding forward** | converts shear into normal forces

**Mid height zip closure** for easy adjustment when client is seated | cushion keeps its shape, enhancing the user's skin protection

**Thigh support**

**PRESSURE INJURY TREATMENT PLAN**

**100% BREATHABLE**

**SOLID SHAPE**

**MACHINE WASHABLE**

**READY FOR USE**

**145 STANDARD SIZES**

- Individual solution for asymmetry
- Minimizes sliding and shear forces
- Optimal adjustment to the body's contours for equal pressure distribution

The Vicair Adjuster 02 consists of 5 individual compartments, designed to contour under loading. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials, coccyx and trochanters, positioning the pelvis and supporting the femurs. For specific needs the filling of each of the 5 compartments can be easily adjusted. Allowing you to change the density of each compartment benefitting users who: (1) need accommodation of a mild to severe fixed pelvic obliquity or correction of a flexible pelvic obliquity; (2) benefit from additional support of the thighs, for instance users with an above knee amputation.

### AVAILABLE COVERS

Vicair 02 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.

### PRODUCT PERFORMANCE

#### Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲▲
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

#### Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲△

#### Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

#### Functionality & Application

Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲▲▲
Maximum user weight	▲▲▲▲△

### TECHNICAL SPECIFICATIONS

Thickness	6 and 10 cm
Average weight	725 grams
Max. user weight	250 kg
Machine washable	60°C max.

### CUSHION SIZES

- Standard (Sizes >53 cm thickness 10cm only)
- On request

Inch	width															
	14	15	16	17	18	19	20	21	22	23	24					
CM	35	38	40	43	45	48	50	53	55	58	60					
14	35															
15	38															
16	40															
17	43															
18	45															
19	48															
20	50															
21	53															
22	55															
23	58															
24	60															

070069-OAD-EN-REV01/2024