



Wheelchair Cushions



	VIGOUR		SPEX CLASSIC		SPEX CLASSIC SUPERHIGH CONTOUR	CONSTRUCTA SHAPE			CONSTRUCTA FLEX			FLEX		XLELLA BARIATRIC		
	STANDARD	HIGH	STANDARD	HIGH		STANDARD	HIGH	SUPERHIGH	STANDARD	HIGH	SUPERHIGH	STANDARD	HIGH	STANDARD	HIGH	
DESCRIPTION	Designed to create the optimum pelvic foundation for the wheelchair user.		Designed to create the optimum pelvic foundation for the wheelchair user.		Creates an optimum pelvic foundation for the user that requires extra pelvic and hip immersion.	A modular cushion which can adapt to accommodate complex postural needs.			A modular cushion which can adapt to accommodate most complex postural needs. With an array of positioning features and enhanced adjustment range.			A truly adaptable and easy approach for even the hardest-to-solve cushion requirements for users with complex postures.		Purpose-designed bariatric cushion range for sizes 22" and above, with a robust design, refined shaping for larger users, and a sturdy foam construction.		
SIZE RANGE (WHEELCHAIR WIDTH" X DEPTH")	8" x 8" to 20" x 22"		8" x 8" to 20" x 22"		8" x 8" to 20" x 22"	14" x 14" to 20" x 22"			14" x 14" to 20" x 22"			12" x 14" to 20" x 20"		22" x 22" to 26" x 28"		
CUSHION HEIGHT (AT CENTRE OF THIGH TROUGH)	8-10": 5cm 11-15": 6.5cm 16-17": 7cm 18-20": 7.5cm		8-10": 5cm 11-15": 6.5cm 16-17": 7cm 18-20": 7.5cm		8-10": 5cm 11-15": 6.5cm 16-17": 7cm 18-20": 7.5cm	14-16": 7cm 18-20": 9cm			14-16": 12cm 18-20": 13cm			14-16": 12cm 17-20": 13cm		9.5cm		
TROUGH HEIGHT (cm)	3cm	5cm	3cm	5cm	Width < 13": 8cm Width > 13": 10cm	3cm	5cm	8cm	3cm	5cm	8cm	3cm	5cm	3cm	5cm	
PRESSURE CARE LEVEL	Medium		Medium		Medium	High			High			Medium		Medium		
ADJUSTABILITY	Non-adjustable		Strategic Positioning Base (foam contouring pads)		Strategic Positioning Base (foam contouring pads)	Strategic Positioning Base (foam contouring pads), adjustable thigh and pelvic supports, adjustable ischial well.			Puzzle joint allowing lateral hip angle and leg length adjustment, adjustable hip angle, adjustable thigh and pelvic supports, adjustable ischial well.			Strategic Positioning Base (foam contouring pads & wedges)		Strategic Positioning Base (foam contouring pads)		
ADJUSTABLE FOR ASYMMETRIC NEEDS	No		Medium adjustability		Medium adjustability	High adjustability, accommodates pelvic obliquities, windsweeping, abduction.			Very high adjustability, accommodates flexion limitations, pelvic obliquities, windsweeping, abduction, leg length discrepancies, internal/external rotation.			High adjustability, accommodates hip flexion asymmetry, hip abduction, hip adduction and internal/external rotation.		Medium adjustability		
STABILITY	Medium	Medium	High	High	High	Medium			Medium			Medium	Medium	High		
CUSTOMISATION OPTIONS	Leg length discrepancy cut-out; Front edge cut-back for leg contractures; Removal of medial thigh support; Custom medial thigh support; Custom lateral thigh support; Pressure-relieving overlay; Honeycomb overlay; Custom sizes. Refer to Spex Script Form and 'Build-Your-Own Cushion' Script Form.					A range of customisations such as leg length discrepancy cut-out, and removable medial thigh support, can be achieved through adjusting the cushion.					Leg length discrepancy cut-out; Front edge cut-back for leg contractures; Removal of medial thigh support; Custom medial thigh support; Custom lateral thigh support; Pressure-relieving overlay; Honeycomb overlay; Custom sizes. Refer to Spex Script Form and 'Build-Your-Own Cushion' Script Form.					
OPTIONS AND UPGRADES	Spare Outer or Incontinence Covers, Rigidiser Base panel, 5° Wedge.					Spare Outer or Incontinence Covers, Rigidiser Base panel, 5° Wedge, 7.5° wedge, Flexy Contour, Extra Pelvic Adjustment Kit, Extra Strategic Positioning Kit.					Spare Outer or Incontinence Covers, Rigidiser Base panel, 5° Wedge.					
WEIGHT LIMIT	WIDTH <		15": 90kg		15": 90kg	15": 90kg			16": 90kg (200 lbs)			16": 90kg (200 lbs)		16": 90kg		N/A
	WIDTH >		15": 160kg		15": 160kg	15": 160kg			18": 136kg (300 lbs)			18": 136kg (300 lbs)			16": 160kg	

Accessories Available:



Universal Positioning Base for use with Gel or Air cushions



Seat Base



Cushion Rigidiser



Spare Covers



Anterior Pelvic Belts

Spex is as diligent as possible in compiling this information. Product availability may differ between regions, please contact your local Product Specialist to confirm availability. Information is correct at time of printing/release: 02/07/24. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Spex. Product specifications are subject to change without notice. © Spex 2024