


# Specification Chart

## CUSHION SIZE CHART



- Red squares are standard sizes.
- Non-standard or custom sizes are available on request.

## PRODUCT DETAILS

Width	Capacity	Product Weight
10"-15"	200lbs (90kg)	(16x16)  2.2kg
16"-20"	300lbs (136kg)	

- Note: This is the weight of the back support only and does not include the hardware. For hardware weight, refer to last page.

## PRODUCT PERFORMANCE

Postural Control				4	
Lateral Stability					5
Postural Asymmetries					5
Changing Needs					5
Pressure Care					5

- Note: This is a guide only and ratings are dependent on correct use of product