

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

**Karma**

KP25.2



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



The KP-25.2 is suitable for medium- and short-distance travel. Lightweight and compact, it features a foldable frame and a removable compact battery pack. Perfect for your daily use, the KP-25.2 is both convenient and portable.

One of the daily difficulties power wheelchair users encounter is manoeuvring in and through tight spaces like doorways and elevators. The KP25.2 is a compact medium-size wheelchair, meaning it is easier to manoeuvre around tight corners and through narrow doorways.



Features and Benefits

- Suitable for mid-short distance journeys.
- Suitable for indoor and outdoor use.
- **No tools required:** The lightweight, foldable frame and the quick release battery design make transport easy.
- **Height adjustable armrests:** According to the user's body requirements, the armrest height can be adjusted to four different heights, allowing for a more comfortable fit.
- **Side transfer:** Side transfer is made easier with the flip-up armrests and the detachable swing-away footrests.
- **Crash Test Approved:** The KP-25.2 is approved for transportation in accordance with ISO 7176-19. This means it is suitable to be used as a seat in a motor vehicle for transportation.

Specifications

Total Length	985mm – 1115mm
Total Width	590mm   640mm   690mm
Total Height	960mm – 1035mm
Seat Width	410mm   460mm   510mm
Seat Depth	410mm   460mm   460mm
Seat Height	550mm
Armrest to Seat Height	220mm – 280mm
Backrest Height	450mm – 520mm
Length   Folded	755mm
Width   Folded	380mm   395mm   395mm
Height   Folded	910mm
Weight without Batteries	35kg   35kg   36kg
Weight with Batteries	67kg   67.6kg   67.6kg
SWL	135kg
Warranty	Wheelchair: 2 years   Batteries: 6 months

Options

Swing away control mount | Elevating Legrests | Stump Support | Folding headrest | Oxygen Holder | USB Charger | Tension Adjustable Backrest



No Tools required



Height Adjustable Armrests



Swing-away Footrests



Notes

Lined area for notes.

# rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider