

A better
wheelchair.
A better fit.



Karma KP-25.2

CODE KP25.2



Frame	Aluminium Folding
Max speed	10km/h
Turning radius	90cm
Travel range	30km
Charger	8 Amp
Batteries	40 ah x 2
Footrest	Swing away removable Height adjustable
Armrest	Swing back Height adjustable Ergonomic pads
Wheels	12" x 2 1/4" mag Pneumatic tyres
Castors	Solid 8" x 2"
Brakes	Electronic brake Push to lock manual brake
Seat	Fire retardant upholstery
Back	Fire retardant upholstery Height adjustable Folding
Drive	Left hand and right hand
Warranty	2 Years - Frame



The KP-25.2 is a foldable electric wheelchair that's comfortable, versatile and adjustable. It's designed to provide easy access and removal of the batteries, making transport even more convenient.

Its flip-back and height adjustable armrests, swing away legs rest assure easy transfers and it is a reliable and easy to use powered wheelchair.

OPTIONS⁺

Swing away control mount | Elevating Legrests | Stump Support | Folding headrest | Oxygen Holder | USB Charger | Tension Adjustable Backrest

SPECIFICATIONS

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight*
KP25.2 - 16"x16"	410mm	410mm	500mm/540mm	480mm	135kg	35kg
KP25.2 - 18"x18"	460mm	460mm	500mm/540mm	480mm	135kg	35kg
KP25.2 - 20"x18"	460mm	460mm	500mm/540mm	480mm	135kg	35kg

*without batteries

rehabhire & sales

310 Lorimer Street Port Melbourne Victoria 3207

t: 1300 000 030 | e: contact@rehabhire.com.au | rehabhire.com.au