
ARJOHUNTLEIGH

GETINGE GROUP

EASYTRACK FS

FEATURING THE V3
PORTABLE LIFT



A true portable lift and transfer system for the home that's easy to transport and quick to assemble.


The Easytrack FS is a freestanding, versatile solution for all types of point-to-point transfers - wheelchair to bed, bed to chair, even in the bathroom, anywhere there is a lift and transfer need.

The *Easytrack FS* has a safe working load of 440 lbs (200 kg). That's exceptional for a portable product making it suitable for use in most any home.

This equipment was designed to allow those with a mobility impairment to continue to live at home and enjoy many of their normal activities.

...with people in mind

EASY ASSEMBLY



The light-weight aluminium bases feature adjustable levels to ensure that the unit remains perfectly stable and secure when assembled.



The posts have been designed to be self-aligning. Components only fit together one way to prevent any assembly errors and ensure the safety of the system.




The rail can be extended up to a width of 115" (3 meters) which is more than enough to span a king size bed.



The telescoping posts include height indicators that allow you to adjust both posts for a perfectly level top rail. A special one-way locking mechanism ensures that the post cannot retract on its own.





- For most transfers
- Head-rest incorporated
- Suitable for most individuals
- Excellent for bathing



- · For general transfers
- Perfectly suited for stocky or bariatric individuals with strong hips or large thighs



- Ideal for transferring to the toilet
- Permits easy access to change and dress individuals

Product			
Easytrack FS Gantry		Part I	Number: 9720010 (*Lift sold separately)
Slings (*Sold separately)	Small	Medium	Large
Weight	45-100 lbs (20-45 kg)	100-200 lbs (45-90 kg)	200-600 lbs (90-272 kg)
Height	4'-5'11" (120-180 cm)	5'-5'11" (152-180 cm)	6'+ (181 cm +)
Hammock Sling	THA-S	THA-M	THA-L
Hammock 6 Sling	THA6-S	THA6-M	THA6-L
Quick-fit Sling	TIR-S	TIR-M	TIR-L
Hygienic Sling	THY-S	THY-M	THY-L
Optional protective carrying bags available - kit #002.17010 (2 post bags, 1 base bag and 1 rail bag) or sold individually			



Use only ArjoHuntleigh designed parts, i.e. parts specifically designed for the purpose, on equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development, we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2011