

ICARE EXTRA LONG IC333

The same functions, features and durability as the standard size IC333.

Overall Bed Length: 2450mm
Extension Bolster Length: 300mm
Please note: The extension and bolster does not function with the knee break.



KEY FUNCTIONS




Head Lift



Knee Break



Height Adjustable
22-66cm



Trendelenburg and
Reverse Trendelenburg

ICARE EXTRA LONG BOLSTERS


Mattress Bolsters are finished in the same fabric as the ActiveX™ Mattresses.
Bolsters available in: 15cm or 20cm heights

IC15 MATTRESS BOLSTER

The 15cm option is compatible with the IC15, M1 and M2 mattresses.

IC20 MATTRESS BOLSTER

The 20cm option is compatible with the IC20, IC30, VersaFlo® and M3 mattresses.



BED CODES & SIZES

LONG SINGLE EXTRA LONG

Code: IC333ELLS
Mattress Surface: 2030 (+300) x 900
Base Size: 2450 x 960

15cm Bolster
Code: MBLS15

20cm Bolster
Code: MBLS20

KING SINGLE EXTRA LONG

Code: IC333ELKS
Mattress Surface: 2030 (+300) x 1070
Base Size: 2450 x 1120

15cm Bolster
Code: MBKS15

20cm Bolster
Code: MBKS20

LONG DOUBLE EXTRA LONG

Code: IC333ELD
Mattress Surface: 2030 (+300) x 1350
Base Size: 2450 x 1400
• One hand control

15cm Bolster
Code: MBLD15

20cm Bolster
Code: MBLD20

ONE PIECE QUEEN EXTRA LONG

Code: IC333ELQ
Mattress Surface: 2030 (+300) x 1520
Base Size: 2450 x 1580
• One hand control

15cm Bolster
Code: MBQ15

20cm Bolster
Code: MBQ20



ASSESSORS NOTE:

Extra Long flat sheets are 300mm longer than Icare Adjustable Bed Sheet Sets. Set includes separate fitted sheet for bolster. Available in white only.