

ICARE IC333

The IC333 is a four function bed in a range of sizes, colours and accessories. This innovative design combines aesthetic appeal with all the essential features. With its ultra-low height, four-functions, and ample under-bed clearance, this bed offers the perfect blend of aesthetics and functionality.



KEY FUNCTIONS




Head Lift



Knee Break



Height Adjustable
22-66cm




Trendelenburg and
Reverse Trendelenburg


COLOUR CHOICES



ONYX



STONE



SMOKE




MIST



ROYAL



NOIR



SILVER



VINYL EBONY
(wipeable)



VINYL DOVE
(wipeable)



SAND
(wipeable)



BISON
(wipeable)




SLATE
(wipeable)

CHOICE OF HEADBOARD/FOOTBOARD DESIGN



3 BUTTON WITH CURVE



STANDARD



KARLSTEAD




AMSTERDAM




COLOGNE

BED CODES & SIZES




LONG SINGLE

Code: IC333LS




ONE PIECE QUEEN

Code: IC333Q




KING SINGLE

Code: IC333KS




SPLIT QUEEN

Code: IC333SQ



LONG DOUBLE

Code: IC333D



SPLIT KING

Code: IC333LS x2

NOTE: Always remember that TWO Long Single beds together make a King size bed.

KEY FEATURES



Extra Long IC333 Option Available



200kg SWL on Long Single & King Single



Unlimited Under Bed Clearance



Available with Headboards and Footboards



External Rechargeable Battery Back-up



Lockable Heavy Duty Castors



300kg SWL on Double and Queen



End and Corner Mattress Retainer Bars



10 Year Warranty



Backlit Wipeable Hand Piece



No gap between mattresses on Queen & King Beds



Strong Breathable Mattress Platform



Ultra Low Design Reducing Risk of Falls



Optional Three Sided Surround for Wheelchair Transfers



Optional IV Pole Attachment Available

COMPATIBLE ACCESSORIES AVAILABLE



High Side Rail

Low Side Rail

Bed Stick

U-Assist Rail

Over Bed Pole

Full Length Side Rails

Headboard/Footboard

Slide-Out Side Rail

Over Bed Tables

Extra Long Bolsters

Bedside Pouch

Bedside Folding Shelf

Bedside Tables

Blanket Cradle