



## CONFIGURA COMFORT ASSESSMENT FORM

Prescriber's Name: \_\_\_\_\_ Client's Name: \_\_\_\_\_

Number: \_\_\_\_\_ Address: \_\_\_\_\_

Email: \_\_\_\_\_ Number: \_\_\_\_\_






Funding Body: \_\_\_\_\_

AREAS HIGHLIGHTED IN GREY ARE STANDARD SPECIFICATIONS



### Size

Chair Size	Part Number	Seat Height	Seat Width	Seat Depth	Weight Capacity
<input type="checkbox"/> SMALL	<input type="checkbox"/> <b>CR5404</b> - DURATEK <input type="checkbox"/> <b>CR5422</b> - Beige Vinyl <input type="checkbox"/> <b>CR5443</b> - Black Vinyl	<input type="checkbox"/> 375 mm (15") #5	412 mm (16.5") #6	<input type="checkbox"/> 350 mm (14") #1	Seat: 160 kg Legrest: 60 kg
		<input type="checkbox"/> 400 mm (16")		<input type="checkbox"/> 400 mm (16") #2	
		<input type="checkbox"/> 425 mm (17") #3	450 mm (18")	<input type="checkbox"/> 450 mm (18")	
		<input type="checkbox"/> 450 mm (18") #4		<input type="checkbox"/> 500 mm (20")	
<input type="checkbox"/> MEDIUM	<input type="checkbox"/> <b>CR5424</b> - Beige Vinyl <input type="checkbox"/> <b>CR5444</b> - Black Vinyl	<input type="checkbox"/> 425 mm (17") #5	462 mm (18.5") #6	<input type="checkbox"/> 350 mm (14") #1	Seat: 160 kg Legrest: 60 kg
		<input type="checkbox"/> 450 mm (18")		<input type="checkbox"/> 400 mm (16") #2	
		<input type="checkbox"/> 475 mm (19") #3	500 mm (20")	<input type="checkbox"/> 450 mm (18")	
		<input type="checkbox"/> 500 mm (20") #4		<input type="checkbox"/> 500 mm (20")	
<input type="checkbox"/> LARGE	<input type="checkbox"/> <b>CR5426</b> - Beige Vinyl <input type="checkbox"/> <b>CR5445</b> - Black Vinyl	<input type="checkbox"/> 425 mm (17") #5	512 mm (20.5") #6	<input type="checkbox"/> 350 mm (14") #1	Seat: 160 kg Legrest: 60 kg
		<input type="checkbox"/> 450 mm (18")		<input type="checkbox"/> 400 mm (16") #2	
		<input type="checkbox"/> 475 mm (19") #3	550 mm (22")	<input type="checkbox"/> 450 mm (18")	
		<input type="checkbox"/> 500 mm (20") #4		<input type="checkbox"/> 500 mm (20")	
<input type="checkbox"/> EXTRA LARGE	<input type="checkbox"/> <b>CR5446</b> - Black Vinyl	<input type="checkbox"/> 375 mm (15") #5	562 mm (22.5") #6	<input type="checkbox"/> 350 mm (14") #1	Seat: 255 kg Legrest: 65 kg
		<input type="checkbox"/> 400 mm (16")		<input type="checkbox"/> 400 mm (16") #2	
		<input type="checkbox"/> 425 mm (17") #3	600 mm (24")	<input type="checkbox"/> 450 mm (18")	
		<input type="checkbox"/> 450 mm (18") #4		<input type="checkbox"/> 500 mm (20")	

#1 Depth Adjustment	#2 Depth Adjustment	#3 Height Adjustment	#4 Height Adjustment	#5 Height Reduction
<input type="checkbox"/> 2 x 2" Kits (100mm) <b>CA2417</b>	<input type="checkbox"/> 2" Kit (50mm) <b>CA2416</b>	<input type="checkbox"/> 1" Raising Kit (25mm) <b>CA2421</b>	<input type="checkbox"/> 2" Raising Kit (50mm) <b>CA2422</b>	<input type="checkbox"/> 1" Lowering Kit (18mm) <b>CA2423</b>
				



Legrest Setting	
Standard - 1" (-25mm)	<input type="checkbox"/>
Standard	<input type="checkbox"/>
Standard + 1" (25mm)	<input type="checkbox"/>
Standard + 2" (50mm)	<input type="checkbox"/>

Armrests	
Overhang on outside (std)	<input type="checkbox"/>
#6 Overhang on inside - reduces width between the armrests by 1.5" (38mm)	<input type="checkbox"/>



Extra Backrest Pillows			
Beige SC2308*	<input type="checkbox"/>	Headrest Pillow	<input type="checkbox"/>
Black SC2565*	<input type="checkbox"/>	Lumbar Pillow	<input type="checkbox"/>

Supports a taller torso by providing 250mm extra height. S, M, & L models only  
\*Full Set of 3 Pillows



## Pressure







Cushion cutout size 450 mm x 450 mm x 75 mm (18"x18"x3")





## Posture



## Cognitive

Waterfall Backrest (std)	Profiled Headrest	Cocoon Backrest **
<input type="checkbox"/> Beige 	<input type="checkbox"/> CA2419 Beige 	<input type="checkbox"/> CA2713 Beige 
<input type="checkbox"/> Black 	<input type="checkbox"/> CA2530 Black 	<input type="checkbox"/> CA2711 Black 



Lateral Support Backrest	Adjustable Lateral Support **
Beige <input type="checkbox"/> CA2406 - SMALL <input type="checkbox"/> CA2407 - MEDIUM <input type="checkbox"/> CA2408 - LARGE 	<input type="checkbox"/> CA2410 - SMALL <input type="checkbox"/> CA2411 - MEDIUM <input type="checkbox"/> CA2412 - LARGE 
Black <input type="checkbox"/> CA2540 - SMALL <input type="checkbox"/> CA2541 - MEDIUM <input type="checkbox"/> CA2542 - LARGE, XL 	<b>Lateral (Pairs)</b> <input type="checkbox"/> CA2413 



\*\* Not compatible with XL size





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



**ENABLE LIFECARE PTY LTD**

**T 1300 370 370**

**[enablelifecare.com.au](http://enablelifecare.com.au)**