

# Metro - Attendant Propelled

Nimble | Compact | Portable



## FEATURES & BENEFITS

- **Compact Folding:** Legrests remain attached to the chair to ensure essential components are not misplaced during transit.
- **Generous Cushioning:** Plush padded mesh seat and backrest provide enhanced seating comfort.
- **Ergonomic Design:** Easy and personalised legrest height adjustment.
- **Flip-Up Armrests:** Simplifies transfers whilst improving access to tables and work surfaces.
- **Carer Confidence:** Equipped with attendant handbrakes for enhanced safety.
- **Essentials Pouch:** Stow small items neatly in the integrated backrest pouch.
- **Oversize Wheels:** 14" rear transit wheels improves ride comfort on uneven terrain.

	OVERALL LENGTH (mm)	OVERALL WIDTH (mm)	OVERALL HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	SEAT HEIGHT (mm)	BACKREST HEIGHT (mm)	ARMREST HEIGHT (mm)	PRODUCT WEIGHT (kg)	REAR WHEEL (inch)	FOLDED DIMENSIONS (mm)
<b>Metro</b> Attendant Propelled MWS442350	930	600	910	400	450	450	470	220	12	14"	680 x 220 x 670