

BED SIDE RAILS

Overview

Explore our bed rail range with clear features and practical applications to help you choose the right support for safer movement, easier transfers, and greater peace of mind at home.

Bed Stick

Features

- Safety return handle
- Mount anywhere on the bed sides

Applications

- Users can use it to pull up, reposition, or sit up.
- Helps users transfer in and out of bed safely



U-Assist

Features

- Slim universal design
- Multiple grab points
- Mount anywhere on the bed sides

Applications

- Multiple handholds for transfers, repositioning, or sitting up.
- Lets users choose the most comfortable or reachable rail for easier movement.




Full Length

Features

- Four rail design
- Caregiver operated
- Left and right mounting options

Applications

- Prevents accidental roll-out, helping users feel secure during sleep or rest.
- Easily folded by caregivers for access or transfers.



Extendable

Features

- Two position extendable
- User operated

Applications

- User can easily retract to allow safe bed entry and exit.
- Provides a secure barrier to reduce the risk of accidental roll-out when extended.



Low

Features

- Hybrid grab / safety rail
- Mount anywhere on bed sides

Applications

- Multiple horizontal rails support transfers and sitting up safely.
- Provides a secure barrier to reduce the risk of roll-out when extended.



High

Features

- Four rail hybrid design
- Mount anywhere along bed sides

Applications

- Ideal for getting in and out of bed.
- Small enough to use as a grab rail, yet large enough to aid positioning and prevent roll-out.



* All rails are designed to IEC60601-2-52 standards.



BED SIDE RAILS

Bed Stick

Our simple grab rail for transfers and repositioning. Slip-on design with safety return hook, adjustable anywhere along the bed.

Height above mattress platform: 520mm

Space inside hook: 103mm

Max user weight: 200kg

Max Static Load: 250N


ARTG Identifier: 259891

SKU: ACBS

Colour Options:



ACBS



Mesh Rail Cover

Breathable padded bed stick cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof + Vant Mesh**



Padded Cover

Waterproof padded bed stick cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof fabric**



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

**Material: Woven fabric
or Vinyl**

Mist ACSRBBI	Noir ACSRBNO	Dove ACSRBVD	Ebony ACSRBVE
Silver ACSRBSI	Stone ACSRBST	Smoke ACSRBSM	Onyx ACSRBON
Bison ACSRBBI	Sand ACSRBSA	Slate ACSRBSL	

*Must be used over padded or mesh rail cover for snug fit.



BED SIDE RAILS

U-Assist

Our slimmest grab rail with multiple grip points for bed transfers. Slip-on design, compatible with Icare homecare beds, mounts anywhere along the side.

Length: 150mm

Height above mattress platform: 530mm

Rail spacing: 100mm

Max user weight: 200kg

Max Static Load: 250N

ARTG Identifier: 259891

SKU: ACUSR

Colour Options:



Charcoal

ACUSR



Mesh Rail Cover

Breathable padded u-assist cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof + Vant Mesh**



ACURCV



Padded Cover

Waterproof padded u-assist cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof fabric**



ACSRCB



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

**Material: Woven fabric
or Vinyl**

Mist ACSRSBMI	Noir ACSRSBNO	Dove ACSRSBVD	Ebony ACSRSBVE
Silver ACSRSBSI	Stone ACSRSBST	Smoke ACSRSBSM	Onyx ACSRSBON
Bison ACSRSBBI	Sand ACSRSBSA	Slate ACSRSBSL	

*Must be used over padded or mesh rail cover for snug fit.



icaremedicalgroup.com

BED SIDE RAILS

Full Length

Our full-length fold-down grab rail for left or right side. Carer-release design with four-bar safety that folds toward foot.

Length: 1385mm

Height above mattress platform: 495mm

Rail spacing: 80mm

Max user weight: 200kg

Max Static Load: 250N

ARTG Identifier: 259891

SKU: ACFSR

Colour Options:



Charcoal

ACFSRR

Right Side

ACFSRL

Left Side

ACFSRP

L+R Pair



Mesh Rail Cover

Breathable padded full length cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof + Vant Mesh**



ACSRCFLV



Padded Cover

Waterproof padded full length cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof fabric**



ACSRCFL



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

**Material: Woven fabric
or Vinyl**

Mist	Noir	Dove	Ebony
ACSRNFLMI	ACSRNFLNO	ACSRNFLVD	ACSRNFLVE
Silver	Stone	Smoke	Onyx
ACSRNFLSI	ACSRNFLST	ACSRNFLSM	ACSRNFLON
Bison	Sand	Slate	
ACSRNFLBI	ACSRNFLSA	ACSRNFLSL	

*Must be used over padded or mesh rail cover for snug fit.



icaremedicalgroup.com

BED SIDE RAILS

Extendable

Our extendable side rail for either side. Slip-on design with four-bar safety, locks in standard or extended length, and positions anywhere along the bed.

Length: 690mm - 950mm

Height above mattress platform: 450mm

Rail spacing: 100mm

Max user weight: 200kg

Max Static Load: 250N

ARTG Identifier: 259891

SKU: ACESR

Colour Options:



Mesh Rail Cover

Breathable padded extendable side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

Material: Zerotec Wipeable Waterproof + Vant Mesh



*Designed to fit snug when fully extended.



Padded Cover

Waterproof padded extendable side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

Material: Zerotec Wipeable Waterproof fabric



*Designed to fit snug when fully extended.



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

Material: Woven fabric or Vinyl

Mist ACSRSEMI	Noir ACSRSENO	Dove ACSRSEVD	Ebony ACSRSEVE
Silver ACSRSESI	Stone ACSRSEST	Smoke ACSRSESM	Onyx ACSRSEON
Bison ACSRSEBI	Sand ACSRSESA	Slate ACSRSESL	

*Must be used over padded or mesh rail cover for snug fit.



icaremedicalgroup.com

BED SIDE RAILS

Low Side

Our support rail for stability and fall prevention. Slip-on design with three-bar safety, positions anywhere along the bed.

Length: 690mm

Height above mattress platform: 422mm

Rail spacing: 115mm


Max user weight: 200kg

Max Static Load: 250N

ARTG Identifier: 259891


SKU: ACLSR

Colour Options:



Charcoal

ACLSR



Mesh Rail Cover

Breathable padded low side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof + Vant Mesh**



ACSRCLV



Padded Cover

Waterproof padded low side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof fabric**

Blue

ACSRCL



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

**Material: Woven fabric
or Vinyl**

Mist ACSRSLMI	Noir ACSRSLNO	Dove ACSRSLVD	Ebony ACSRSLVE
Silver ACSRSLSI	Stone ACSRSLST	Smoke ACSRSLSM	Onyx ACSRSLON
Bison ACSRSLBI	Sand ACSRSLSA	Slate ACSRSLSL	

*Must be used over padded or mesh rail cover for snug fit.



icaremedicalgroup.com

BED SIDE RAILS

High Side

Our support rail for safe bed entry and exit. Slip-on design with four-bar safety, fits Icare healthcare beds, positions anywhere along the side.

Length: 350mm

Height above mattress platform: 582mm

Rail spacing: 119mm

Max user weight: 200kg

Max Static Load: 250N

ARTG Identifier: 259891

SKU: ACHSR

Colour Options:



Mesh Rail Cover

Breathable padded high side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof + Vant Mesh**



Padded Cover

Waterproof padded high side rail cover to protect against hard impacts and prevent limbs from slipping between the rails.

**Material: Zerotec Wipeable
Waterproof fabric**



Fabric Slip Cover

Fabric slip covers fit over the mesh or padded rail covers and can match the bed surround upholstery.

**Material: Woven fabric
or Vinyl**

Mist ACSRSHMI	Noir ACSRSHNO	Dove ACSRSHVD	Ebony ACSRSHVE
Silver ACSRSHSI	Stone ACSRSHST	Smoke ACSRSHSM	Onyx ACSRSHON
Bison ACSRSHBI	Sand ACSRSHSA	Slate ACSRSHSL	

*Must be used over padded or mesh rail cover for snug fit.

