
ARJOHUNTLEIGH

MAXI MOVE

A MOBILE LIFTER FOR BUSY CARE ENVIRONMENTS

GETINGE GROUP



Maxi Move[™], the most complete mobile lifter system available, is a proven success in care environments around the world. ArjoHuntleigh has now enhanced the *Maxi Move* to make it an even greater asset for care settings of all types.

The fundamental aim of *Maxi Move* is to be "best on basics". Using *Maxi Move* you gain a safer, more efficient solution for the basic tasks of resident handling – lifting and repositioning.

The system offers two unique features – SVS (Stable Vertical System) and Powered DPS (Dynamic Positioning System), which together provide the safest, most comfortable basis for complete transfers – from the initial lift through to effortless fine-tuning of the resident's posture.

This battery-powered passive sling lifter allows a single carer to safely manage a wide range of patient handling routines. With unrivalled flexibility, stability and reliability, *Maxi Move* is the best-equipped lifter for intensive use in busy care environments.












FLEXIBILITY FOR THE DIVERSE CHALLENGES OF EVERYDAY CARE

There is no such thing as a typical patient or a typical care facility. That is why *Maxi Move* has been designed as a complete system with the flexibility to provide patient-specific solutions in any care environment.





An optimum solution for supine lifting of patients in acute care...



...rehabilitation sessions using the unique Walking jacket...



...or for repositioning and turning patients in bed.



A safe and easy way to position patients in bed...



...transfer patients to chairs...



...or lift supine patients on stretchers.




ArjoHuntleigh is a branch of Arjo Ltd Med. AB. Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.
© and TM are trademarks belonging to the ArjoHuntleigh group of companies.
© ArjoHuntleigh, 2010

1 Scale

Built-in electronic scale shows resident's weight on the display of hand control.



The Combi attachment enables easy interchanging between spreader bars and stretchers.



DPS (Dynamic Power System) and Power DPS minimises manual handling and improves precision.



User friendly hand control with a colored and illuminated screen.

Accessories

A wide selection of spreader bars

Wide range of special purpose slings and flites

Stretcher solutions

Optional: extended jib,

scale (also available as an upgrade kit)

Arjo Clean (for disinfection)

For tecnical specifications of Maxi Move, please visit www.ArjoHuntleigh.com