



MONET DUAL ACTION

Dual Action

The Monet Dual Action offers independent control of back, leg and headrest functions, giving the user infinite positioning options at their fingertips. The heating function provides warmth and comfort during the colder months whilst Zero Gravity tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.




Features & Benefits

- **Dual Action** - Independent control of leg and recline for infinite positioning.
- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book.
- **Integrated Lift and Recline** functions for easy sit-to-stand.
- **Extending Legrest** supports more of the lower legs for enhanced comfort.
- **Heating Function** within the backrest provides warmth and comfort during the colder months.
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.
- **Backlit Hand Control** with integrated USB charging point.
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day.

5 YEAR WARRANTY
FRAME

2 YEAR WARRANTY
COMPONENTS*

S.W.L
159KG

MONET DUAL ACTION CHP227540	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	BACKREST HEIGHT (mm)	OVERALL HEIGHT (mm)	OVERALL WIDTH (mm)	OVERALL DEPTH (mm)	DISTANCE FROM WALL (mm)
								
	495	535	560	570	1065	890	1065	405

* Components: Handset, upholstery and electronics.