



WALKER RANGE

Premium Mobility Solutions



VOGUE CARBON FIBRE SEAT WALKER

One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.

CARBON FIBRE FRAME: Extremely lightweight to maximise portability

COLOURS:

- Small - Sapphire
- Medium - Black
- Tall - Emerald



MODELS

WAF705340	Vogue Carbon Fibre - Small - Sapphire
WAF705350	Vogue Carbon Fibre - Medium - Black
WAF705370	Vogue Carbon Fibre - Tall - Emerald



LESS THAN 5KG*

VOGUE LIGHTWEIGHT SEAT WALKER

This modern seat walker has a clean and crisp frame design, making it stand out when going for a walk or browsing the shops!

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:

- Red
- Flower
- Champagne
- Zebra
- Leopard

MODELS

WAF705310RE	Vogue Lightweight - Red
WAF705310CH	Vogue Lightweight - Champagne
WAF705310LE	Vogue Lightweight 2 - Leopard
WAF705310FL	Vogue Lightweight 2 - Flower
WAF705310ZB	Vogue Lightweight 2 - Zebra



Discreet Cable Management

WEIGHS JUST OVER 5KG



FEATURES & BENEFITS SHARED BY THE VOGUE CARBON FIBRE ...




Ergonomic Design

Offers personalised height adjustment which promotes correct posture



Compact

Turning circle for improved mobility in tight spaces



Simple

Folding mechanism locks for easy transportation and storage



Handbrakes

Lock for user safety and stability


	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
Small											
	655mm	600mm	715-830mm	240mm	460mm	500mm	230mm	203mm	4.7kg	5.9kg	5kg
Medium											
	26"	23.6"	28"-32"	9.4"	18"	20"	9"	8"	10.4lbs	13lbs	11lbs
Tall											
	655mm	615mm	725-940mm	240mm	460mm	550mm	230mm	203mm	4.95kg	6.2kg	5kg
Medium											
	26"	24"	28"-37"	9.4"	18"	22"	9"	8"	11lbs	13.7lbs	11lbs
Tall											
	655mm	600mm	780-1030mm	240mm	460mm	620mm	230mm	203mm	5.2kg	6.35kg	5kg
Tall											
	26"	23.6"	30"-40"	9.4"	18"	25"	9"	8"	11.5lbs	14lbs	11lbs

... AND LIGHTWEIGHT ALUMINIUM SEAT WALKERS



Storage

Removable bag carries up to 5kg of essentials



Backrest

Deep seat and wide backrest for improved comfort and support



Bonus

Walking cane holder

	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
Small											
	655mm	620mm	785-920mm	240mm	465mm	540mm	220mm	203mm	5.3kg	6.45kg	5kg
Medium											
	26"	24.5"	31"-36"	9.5"	18"	21"	8.6"	8"	11.7lbs	14.2lbs	11lbs



AS / NZS Approved


Vogue is registered in the ARTG, AUST R ARTG 271255

* Carbon Fibre Walker - Medium | ** Includes accessories, bag and backstrap.

VOGUE ADVENTURE SEAT WALKER

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

COLOURS:



Wider Seat supports up to 200kg



Large Multi-Terrain Wheels



Folds to 245mm



MODELS

WAF705450 | Vogue Adventure - Black

FEATURES & BENEFITS

Ergonomic Design
Offers personalised height adjustment which promotes correct posture

Handbrakes
Lock for user safety and stability

Simple
Folding mechanism locks for easy transportation and storage

Compact
Turning circle for improved mobility in tight spaces

Storage
Removable bag carries up to 5kg of essentials

Bonus
Walking cane holder

SWL 200kg 31.5st

1 YEAR WARRANTY

AS / NZS Approved

Vogue is registered in the ARTG, AUST R ARTG 271255


Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size (Front / Rear)	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
770mm	695mm	900-1035mm	230mm	540mm	530mm	245mm	279mm / 225mm	8.8kg	10kg	5kg
30.3"	27.4"	35.4"-40.7"	9.1"	21.3"	20.9"	9.6"	11" / 8.8"	19.4lbs	22lbs	11lbs

VOGUE FOREARM WALKER



For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:



White / Black



Padded Forearm Supports



Folds to 225mm



MODELS

WAF705400 | Vogue Forearm - White/Black

FEATURES & BENEFITS

Ergonomic Design
Offers personalised height adjustment which promotes correct posture

Handbrakes
Lock for user safety and stability

Bonus
Walking cane holder

Simple
Folding mechanism locks for easy transportation and storage

Storage
Removable bag carries up to 5kg of essentials

SWL 150kg 23.6st

1 YEAR WARRANTY

AS / NZS Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
660mm	620mm	1020-1270	240mm	460mm	620mm	225mm	203mm	9.1kg	10kg	5kg
26"	24.4"	40.2"-50"	9.4"	18.1"	24.4"	8.8"	8"	20.1lbs	22lbs	11lbs

** Includes accessories, bag and backstrap.

VOGUE INDOOR WALKER

Ideal for use within the home, the Indoor Walker is perfect for portable dining and in home mobility assistance.

COLOURS:



Red / Black



Champagne / Black

MODELS

WAF705100RE	Vogue Indoor - Red/Black
WAF705100CH	Vogue Indoor - Champagne/Black



Folds to 270mm

FEATURES & BENEFITS



Ergonomic Design
Offers personalised height adjustment which promotes correct posture



Handbrakes
Lock for user safety and stability



Simple
Folding mechanism locks for easy transportation and storage



Portable Dining
Safely collect and transport meals or hot drinks



Compact
Turning circle for improved mobility in tight spaces



Storage
Removable bag carries up to 5kg of essentials



Bonus
Walking cane holder



Turn Around.
Slide Over Knee.
Lock & Enjoy.

SWL
110kg
17.3st

1 YEAR WARRANTY

✓ AS / NZS Approved

✓ Vogue is registered in the ARTG, AUSTRALIA ARTG 271255

Overall Length	Overall Width	Overall Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL	Tray SWL
650mm	540mm	860-970mm	270mm	152mm	4.4kg	5.1kg	5kg	5kg
25.6"	21.3"	33.9"-38.2"	10.7"	6"	9.7lbs	11.2lbs	11lbs	11lbs

** Includes accessories, bag and tray.

VOGUE WALKER ACCESSORIES



Oxygen Bottle Holder Bag

WAA695600



Night Light - LED Lamp

WAA695630



IV Pole Lamp

WAA695651



Single Hand Brake System* - Compatible with Aspire Vogue Carbon Walkers

WAA695640

*Single Hand Brake System are not compatible with the Vogue Indoor Walker



ASSOCIATED PRODUCTS



MOBILITY SCOOTERS

Wide range of Aspire Scooters for all user needs.

Shown Aspire HS268 Manual Folding Mobility Scooter SCT671135



COMFIMOTION COLLECTION

Range of adjustable beds and accessories for both clinical needs and comfort at home.


Shown Aspire ComfiMotion Activ Care Bed in Portebello Grey - Long Single - BEB046740



LIFT RECLINE CHAIRS

Luxury Aspire recliner range has four, dual and single motor functions.

Shown Aspire PostureFit Configurable Lift Recline Chair
Sizes Available Medium & Tall - Medium/Large



OVERBED TABLES

Lightweight mobile table positions easily over the lifestyle community beds.

Shown Over Bed Table with Gunmetal Frame & Sublime Teak Laminate Top BEA006705