



# WALKER RANGE

Aspire for... Comfort, stability, support



Aspire Walkers provide assistance, support, and relief to people who want to remain actively mobile and independent — giving confidence to users of all shapes, sizes, and mobility levels.

## COMMON FEATURES & BENEFITS

✓ DVA Contracted

✓ AS / NZS Approved

1 YEAR LIMITED WARRANTY




### Rustproof Aluminium Frames

For durability and stability



### Adjustable Handle Heights

Facilitating optimal sizing for postural positioning



### Latex-free Ergonomically Shaped Grip Handles

Incredibly comfortable for users with arthritis



### Ultra-Soft Hand Brakes

With effortless palm ball lock



### Easy Folding Mechanism

Facilitating compact storage and transportation



### Flip up Padded Seat

For accessing clip-on removable storage bags with discreet pocket compartments



### Curved Padded Backrests

For seated support and comfort



### 360° Swivelling Front Wheels

For easy movability



### Easy Wipe

For clean surfaces



### Tool Free

Assembly and adjustment

## CLASSIC 8 SEAT WALKER

**Handle Height Range**  
820 – 975mm

**Taller Seat Height**  
575mm

**8" Solid Wheels**  
Ease over even the most robust indoor and outdoor terrain



Available colours:

■ WAF705700BL ■ WAF705700RE

✓ Classic 8 is registered in the ARTG, AUST R ARTG 224644

## CLASSIC 6 SEAT WALKER

**Handle Height Range**  
760 – 915mm

**Medium Seat Height**  
525mm

**6" Solid Wheels**  
Ease of mobility over every day indoor and outdoor surfaces



Available colours:

■ WAF709600BL ■ WAF709600RE

✓ Classic 6 is registered in the ARTG, AUST R ARTG 224644

## DELUXE 8 SEAT WALKER

**Deluxe, In-Line Rear Wheels and Brakes**  
Increases space around the user's stepping zone

**Handle Height Range**  
835 – 990mm

**Taller Seat Height**  
590mm

**8" Solid Wheels**  
Ease over even the most robust indoor and outdoor terrain



Available colours:

■ WAF705750BL ■ WAF705750RE

✓ Deluxe is registered in the ARTG, AUST R ARTG 224650

## MINI 6 SEAT WALKER

**Handle Height Range**  
745 – 905mm

**Lower Seat Height**  
485mm

**6" Solid Wheels**  
Ease over every day surfaces indoors and out



Available colours:

■ WAF705800BL ■ WAF705800RE

✓ Mini 6 is registered in the ARTG, AUST R ARTG 224651

## FLEX 6 SEAT WALKER

**Handle Height Range**  
730 – 830mm

**Seat Height Range**  
430 – 520mm

**Adjustable** seat and handle height

**6" Solid Wheels**  
Ease of mobility over every day indoor and outdoor surfaces



Available colour:

■ WAF709660CH

✓ Flex 6 is registered in the ARTG, AUST R ARTG 224644

## FLEX 8 SEAT WALKER

**Handle Height Range**  
780 – 970mm

**Seat Height Range**  
475 – 570mm

**Adjustable** seat and handle height

**8" Solid Wheels**  
Ease over even the most robust indoor and outdoor terrain



Available colour:

■ WAF709670BK

✓ Flex 8 is registered in the ARTG, AUST R ARTG 224644

## XL 8 SEAT WALKER

**Handle Height Range**  
850 – 985mm

**Wide Padded Seat**  
With 180kg safe working load for bigger-sized users

**Medium Seat Height**  
560mm

**8" Solid Wheels**  
Ease over even the most robust indoor and outdoor terrain



Available colour:

■ WAF750020SL

✓ XL 8 is registered in the ARTG, AUST R ARTG 224645

## KNEE WALKER

**Mobility Support**  
For an injury below the knee

**Knee Pad**  
Adjustable height

**8" Wheels**  
For outdoor manoeuvrability

**Fold-Down Mechanism**  
For efficient storage



Available colour:

■ WAF712200

✓ Knee Walker is registered in the ARTG, AUST R ARTG 306095

## TRI WHEEL WALKER

**Handle Height Range** 825 to 925mm

**Ideal for use** in confined spaces

**Robust steel** construction

**Folds for easy storage**  
and transportation

**Storage bag**, basket and meal trays

**8" Solid Wheels** provide  
smooth manoeuvrability over  
variable surfaces



Available colour:

 WAF705500



## MEAL TRAY WALKER

**Ideal for transporting meals** around  
the home

**Clip On Plastic Trays** with outer lip  
edge to capture spillage and prevent  
items from slipping off

**Height Adjustable Handles** with  
lockable handbrakes to hold the  
walker in place

**Large 6" Wheels** easier to navigate  
on carpets and door thresholds

**360° Swivelling** front wheels



Available colour:

 DLG272061



## MEAL TROLLEY

**Ideal for transporting meals** around  
the home

**Clip On Plastic Trays** with outer lip  
edge to capture spillage and prevent  
items from slipping off

**3" Castors** with two locking castors  
to hold the trolley in place



Available colour:

 DLG272060



*\*Please note, not all common features and benefits are applicable to these products*

## SPECIFICATIONS

*All measurements shown are in millimetres*

| WALKERS   | OVERALL<br>HEIGHT | OVERALL<br>LENGTH | OVERALL<br>WIDTH | SEAT<br>DEPTH | SEAT<br>WIDTH | SEAT<br>HEIGHT | BACKREST<br>HEIGHT | FOLDED<br>WIDTH | WHEEL<br>SIZE | UNIT<br>WEIGHT | SWL   |
|-----------|-------------------|-------------------|------------------|---------------|---------------|----------------|--------------------|-----------------|---------------|----------------|-------|
| CLASSIC 8 | 820 - 975         | 660               | 630              | 305           | 355           | 575            | 815                | 200             | 8"            | 7kg            | 130kg |
| CLASSIC 6 | 760 - 915         | 625               | 615              | 305           | 355           | 525            | 770                | 200             | 6"            | 6kg            | 130kg |
| DELUXE 8  | 835 - 990         | 765               | 605              | 305           | 355           | 590            | 840                | 245             | 8"            | 7kg            | 130kg |
| FLEX 6    | 730 - 840         | 630               | 550              | 448           | 448           | 436 - 524      | 720                | 220             | 6"            | 6.5            | 136kg |
| FLEX 8    | 777 - 975         | 485               | 589              | 485           | 450           | 483 - 572      | 758                | 293             | 8"            | 7.5            | 136kg |
| MINI 6    | 745 - 905         | 560               | 615              | 305           | 355           | 485            | 690                | 200             | 6"            | 6kg            | 130kg |
| XL 8      | 850 - 985         | 700               | 735              | 340           | 495           | 560            | 850                | 230             | 8"            | 8kg            | 180kg |
| TRI WHEEL | 825 - 925         | 575               | 715              | N/A           | N/A           | N/A            | N/A                | 275             | 8"            | 5.5kg          | 115kg |

| WALKERS     | OVERALL<br>HEIGHT | OVERALL<br>LENGTH | OVERALL<br>WIDTH | KNEE PAD<br>DEPTH | KNEE PAD<br>WIDTH | KNEE PAD<br>HEIGHT | HANDLE<br>HEIGHT | WHEEL<br>SIZE | UNIT<br>WEIGHT | SWL   |
|-------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------------|------------------|---------------|----------------|-------|
| KNEE WALKER | 840 - 970         | xxx               | 430              | 360               | 180               | 480 - 560          | 840 - 970        | 8"            | 10kg           | 145kg |

| MEAL<br>WALKING AIDS | OVERALL<br>HEIGHT | OVERALL<br>LENGTH | OVERALL<br>WIDTH | LOWER TRAY<br>HEIGHT | UPPER TRAY<br>HEIGHT | TRAY<br>WIDTH | WHEEL<br>SIZE | UNIT<br>WEIGHT | SWL   |
|----------------------|-------------------|-------------------|------------------|----------------------|----------------------|---------------|---------------|----------------|-------|
| MEAL TRAY WALKER     | 770 - 920         | 560               | 585              | 415                  | 770 - 920            | 450           | 6"            | 8kg            | 125kg |
| MEAL TROLLEY         | 810 - 940         | 505               | 460              | 355 - 480            | 615 - 740            | 450           | 3"            | 6kg            | N/A   |

## ACCESSORIES



### OXYGEN BOTTLE HOLDER

Attaches to A-Frame  
Seat Walkers  
WAA686905



### SEAT WALKER BASKET

Clamp On Basket To Suit  
Aspire Seat Walkers  
WAA686910 - Deluxe  
WAA686909 - Classic,  
Mini & Push Down  
WAA686908 - Aspire XL



### WALKING STICK HOLDER

Attaches To A-Frame  
Seat Walkers  
Top mount clamp to  
securely hold the walking  
stick in place  
WAA686924