



**MULTiWAY**

Innovation & independence

# ABOUT MULTiWAY

MULTiWAY is an Australian-owned company that prides itself on providing quality products at competitive prices.

We work closely with healthcare clinicians and customers in order to identify the latest technologies and changing needs of the market and our clients.

As specialists in sleep systems, we believe that a bed is an important aspect of your daily care. It is your personal space where quality rest and recovery should be expected, comfortable and rewarding.

The BodyAdjust Sleep System (BASS) is designed to blend into the home environment with ease, yet not compromise on your therapeutic needs. Multiple positioning options will assist you with independent mobility, improved sleep quality and comfort levels.

Available in two versions, Generation II and Generation III




# Generation II



*A quality hi-lo bed with trendelenberg - head and foot elevation.*

## Product Specifications

- Colour  
Midnight black 
- Safe working load  
195kg
- Hi-low lift range  
690mm
- Floor to base (bottom)  
155mm bed in lowest position
- Floor to base (top)  
240mm bed in lowest position
- Castors  
55mm x 8 nylon (4 lockable)
- Trendelenberg tilt angles  
17 degrees
- Warranty  
5 Year - Frame  
2 Year - Electrics
- Maximum back raise angle  
90 degrees
- Maximum leg raise angle  
46 degrees
- Hoist clearance
- Base panels included as standard
- Corded back lit remote control
- Mattress brackets included

**Voice  
ACTIVATION**




\*optional

# Generation II

## Options

- Cordless handset
- Smart device app controller
- Voice control kit
- Massage kit (comes with cordless controller)
- Headboard
- Footboard
- Rubber castor kit (8 lockable)
- 3/4 fold down side rails (pair)
- Side grip handles (pair)

## Sizes



## Product codes

Product Code	Name / Description	Size
BA100	Adjustable Base Only LS	Long Single 915mm x 2030mm
BA101	Adjustable Base Only KS	King Single 1070mm x 2030mm
BA104	Head Board LS	Long Single Headboard 915mm
BA105	Head Board KS	King Single Headboard 1070mm
BA108	3/4 Length Rails	3/4 Length Rails - Pair
BA109	Grip Handle	Grip Handle - Set of 2
BA123	Foot Board LS	Long Single Foot Board 915mm
BA124	Foot Board KS	King Single Foot Board 1070mm
BA400	Voice Activation	Voice Activation (Cordless Handset)
BA411	Massage System	Massage System(Cordless Handset)




# Generation III



*A quality hi-lo bed with trendelenberg - Head and Feet elevation. Included as standard - battery back-up, underbed lighting, usb ports.*

## Product Specifications

- Colour  
Pebble grey 
- Safe working load  
195kg
- Hi-low lift range  
690mm
- Floor to base (bottom)  
155mm bed in lowest position
- Floor to base (top)  
240mm bed in lowest position
- Castors  
55mm x 8 nylon (4 lockable)
- Trendelenberg tilt angles  
17 degrees
- Warranty  
5 Year - Frame  
2 Year - Electrics
- Maximum back raise angle  
90 degrees
- Maximum leg raise angle  
46 degrees
- Hoist clearance
- Base panels included as standard
- Cordless back lit remote control with custom memory setting
- Under bed LED Light
- Lithium Backup Battery
- USB Charge ports (2 x USB-A)
- Mattress brackets included

**Voice  
ACTIVATION**



\*optional

# Generation III

## Options

- Corded handset
- Smart device app controller
- Voice control kit
- Massage kit (comes with cordless controller)
- Headboard
- Footboard
- Rubber castor kit (8 lockable)
- 3/4 fold-down side rails (pair)
- side grip handles (pair)

## Sizes



## Product codes

Product Code	Name / Description	Size
BA100-GEN-III	Adjustable Base Only Long Single	Long Single 915mm x 2030mm
BA101-GEN-III	Adjustable Base Only King Single	King Single 1070mm x 2030mm
BA102-GEN-III	Adjustable Base Only Double	Double 170mm x 2030mm
BA103-GEN-III	Adjustable Base Only Queen	Queen 1525mm x 2030mm
BA104-GEN-III	Head Board Pebble Grey - Single	Head Board Pebble Grey BASS 950mm
BA105-GEN-III	Head Board Pebble Grey - King Single	Head Board Pebble Grey BASS 1070mm
BA106-GEN-III	Head Board Pebble Grey - Double	Head Board Pebble Grey BASS 1370mm
BA107-GEN-III	Head Board Pebble Grey-Queen	Head Board Pebble Grey BASS 1525mm
BA108	3/4 Length Rails	3/4 Length Rails - Pair
BA109	Grip Handle	Grip Handle - Set of 2
BA123-GEN-III	Foot Board Pebble Grey - Single	Foot Board Pebble Grey BASS 950mm
BA124-GEN-III	Foot Board Pebble Grey - King Single	Foot Board Pebble Grey BASS 1070mm
BA125-GEN-III	Foot Board Pebble Grey - Double	Foot Board Pebble Grey BASS 1370mm
BA126-GEN-III	Foot Board Pebble Grey - Queen	Foot Board Pebble Grey BASS 1525mm
BA400-GEN-III	Voice Activation	Voice Activation (Cordless Handset)
BA411-GEN-III	Massage System	Massage System(Cordless Handset)
BA507	Corded Back Lit Handset Controller	Corded Back Lit Handset Controller

# MULTiWAY mattresses

## "Active" inner pocket spring

Product Code	Name / Description	Size
BA110	Active Mattress LS	Long Single 915mm x 2030mm x 230mm
BA111	Active Mattress KS	King Single 1070mm x 2030mm x 230mm
BA112	Active Mattress DBL	Double 1370mm x 2030mm x 230mm
BA113	Active Mattress Q	Queen 1525mm x 2030mm x 230mm

A multi-layered mattress giving you all the benefits of contouring foam, the breathability of cotton and bamboo fibres, and the comfort of inner pocket springs. Ideal for those requiring minimal partner disturbance through the night. The Active also has a dense foam surround to assist with transfer in and out of bed.

- Diamond stitching
- Prevents "bunching"
- Bamboo overlay
- Cooling memory foam gel



# Mattresses

## "Comfort" memory foam - cool gel

Product Code	Name / Description	Size
BA114	Comfort Mattress LS	Long Single 915mm x 2030mm x 230mm
BA115	Comfort Mattress KS	King Single 1070mm x 2030mm x 230mm
BA116	Comfort Mattress DBL	Double 1370mm x 2030mm x 230mm
BA117	Comfort Mattress Q	Queen 1525mm x 2030mm x 230mm

This supple 2-layer combination foam mattress supports and conforms to the contours of your body closely. The bamboo and cool gel upper-layer allows air to circulate and improves immersion to assist with pressure redistribution.

- Bamboo overlay
- Castellated foam
- Inner pocket spring
- Dense foam surround for safer transfers





# Mattress over-lays

MULTiWAY bamboo pillow top overlays provide a little extra softness

Product Code	Name / Description	Size
BA118	Mattress Overlay LS	Long Single 915mm x 2030mm x 230mm
BA119	Comfort Mattress KS	King Single 1070mm x 2030mm x 230mm
BA120	Comfort Mattress DBL	Double 1370mm x 2030mm x 230mm
BA121	Comfort Mattress Q	Queen 1525mm x 2030mm x 230mm

Product Code	Name / Description	Size
BA122	Memory Foam Pillow	One



# Resources

scan to MULTiWAY  
channel



# Resources

scan to MULTiWAY flyer





**MULTiWAY**

**HEAD OFFICE**

**P.** 1300 168 584 (1300 1 MULTI)

**E.** [sales@multiway.com.au](mailto:sales@multiway.com.au) **W.** [multiway.com.au](http://multiway.com.au)