



Ezi-Go DLX


AUSTRALIA

www.meritsaustralia.com


Ezi-Go DLX

Features

- New vibrant fresh colours available
- Dismantle in seconds for easy transportation and storage
- Puncture proof tyres front and rear
- Suitable for indoors and shopping centres
- New stylish comfortable contoured high back seat with height & width adjustable arm rests
- Joystick mounting for left and right hand operation
- Easy to operate quick release trigger which separates the base sections
- Intelligent, regenerative electro magnetic brakes activate automatically when the Joy Stick is released



Dimensions cm



Overall Dimension : A x B x C 90 x 53 x 92cm/35.5 x 21 x 36in

Seat Dimension : D x E x F 46 x 40 x 48cm/18 x 16 x 19in

Seat Height from Ground : G 52-56cm/20.5-22in

Seat Height from Deck : H 42-46cm/16.5-18in

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
25/55	4.5/10	11/24	45/99	114/250

Technical Specs.

Max Speed	6.4kph/4mph
Range	15km/9mi
Turn Radius	57cm/22"
Ground clearance	4.5cm/1.8"
Motor	DC 24V, 90W(rating)/300W(Max)
Controller	PG VR2 50A
Battery	12V/14A/H x 2pcs
Charger	2A off-board
Gradient	6°
Caster	6" x 2" PU tire
Drive Wheel	200 x 75 PU tire
Brake	intelligent, regenerative, electromagnetic brakes

 **merits**
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.