



Maverick 12




Maverick 12

Features

- Mid-wheel drive design offers outstanding manoeuvrability in tight, compact spaces and the superior adjustable suspension design delivers a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension that works
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Swing away joy stick standard
- Height and width adjustable armrests on the new contoured seat. Height adjustable four post seat with reclining back rest
- 204kg weight capacity
- Available in a number of vibrant new colours – Red, Pearl Gold, Charcoal, NV Green, Kite Blue
- Other options available are electric seat lift, rear attendant controls, optional seating & leg rests

Dimensions cm / in



Overall Dimension : A x B x C	90x 69 x 127cm / 35 x 27 x 50in
Seat Dimension : D x E x F	51 x 46 x 51cm / 20 x 18 x 20in
Seat Height from Ground : G	56-58cm / 22-23in
Seat Height from Deck : H	43-48cm / 17-19in

Weight kg / lb

base	battery/ea	seat	total with battery	capacity
51/113	17/37	21/46	106kg /233lb	205kg /450lb



Technical Specs.

Max Speed	8kph / 5mph
Range	32km / 20mi.
Turn Radius	53cm / 21"
Ground clearance	5.6cm / 2.2"
Motor	DC24V, 200W(Rating) / 575W(Max)
Controller	Dynamic Shark 50
Battery	12V / 55-62A/H x 2pcs
Charger	5A off-board
Gradient	8°
Caster	front: 6" x 2" PU tire ; rear: 6" x 2" solid tire
Drive Wheel	12" (4.00-5) pneumatic tire
Brake	intelligent, regenerative, electromagnetic brakes