

RELAX BED

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.



Portobello Grey

*Mattress sold separately



Wireless handset



Back Raise



Leg Raise



Head Tilt



Zero Gravity

FEATURES & BENEFITS

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Press & Hold Safety Handset:** Intuitive, large button interface enhances usability.
- **Box-on-Box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- **USB Charge:** Connect and charge your devices conveniently.
- **Height Adjustable Legs:** Allows the bed height to be adjusted to 320, 370, 380 or 430mm tall. No tools are needed to make the adjustment.

FUNCTIONS

- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- **Leg Raise:** Promotes circulation, relieving tired, swollen and/or painful limbs.
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button.

1 YEAR WARRANTY
HANDSET

2 YEAR WARRANTY
ELECTRONICS

10 YEAR WARRANTY
FRAME



RELAX BED


LONG SINGLE **BEB046900**
KING SINGLE **BEB046910**
LONG DOUBLE **BEB046920**
QUEEN **BEB046930**

OVERALL WIDTH
(mm)



920
1070
1370
1520

OVERALL LENGTH
(mm)



2030

OVERALL HEIGHT
(mm)



Height adjustable to
320, 370, 380 or 430

ACCESSORIES


ALLURE HEADBOARD



- Stylish headboard to fit ComfiMotion bed collection.
- Subtle square padded design.

Long single **BEA830120**
King single **BEA830130**
Long double **BEA830132**
Queen **BEA830134**
King **BEA830136**


CHESTERFIELD HEADBOARD



- Stylish headboard to fit ComfiMotion bed collection.
- Chesterfield padded design.

Long Single **BEA830100**
King Single **BEA830110**
Long Double **BEA830112**
Queen **BEA830114**
King **BEA830116**

MEMORY FOAM MATTRESS



- Integrated side – walls.
- Contoured Immersion.
- Superior breathability.
- Anti- microbial.

Classic – Long Single **BEM469910**
Classic – King Single **BEM469920**
Classic – Long Double **BEM469922**
Classic – Queen **BEM469924**
Incontinence – Long Single **BEM469930**
Incontinence – King Single **BEM469940**
Incontinence – Long Double **BEM469942**
Incontinence – Queen **BEM469944**

POCKET SPRING MATTRESS



- Weight activated support zones that contour to individual body shapes to improve comfort, support and sleep quality.
- Integrated side – walls.
- Luxuriously soft pillowtop layer.

Long Single **BEM469800**
King Single **BEM469810**
Long Double **BEM469820**
Queen **BEM469830**


CROOK HANDLE BED STICK



BEA830200

- Provides multiple grip points to assist repositioning in bed, and safe entry and exit.
- Designed to be used in conjunction with ComfiMotion bed collection.

BED RAIL



BEA830210

- Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.
- Designed to be used in conjunction with ComfiMotion bed collection.

BREEZE PILLOW



BEP050100

- Surface Cooling Technology helps regulate your body temperature for a more comfortable sleep experience.
- Medium Density, Immersive Memory Foam.

PLUSH PILLOW



BEP050110

- Plush Fibrefill sleep surface and Memory Foam core provides ergonomic support and comfort.
- Medium density.