





ICARE IC777

CARING FOR YOUR COMFORT



ICARE IC777 PREMIUM HOMECARE BED



VINYL

Heavy commercial grade jacquard, vinyl fabric.



This four function bed, features dual lift motors totaling 12,000 newtons lifting power, unlimited under-bed clearance and eight heavy duty castors.



62° Head Lift



42° Knee Break



Height Adjustable Min 22 - Max 66cm



Trendelenburg and Reverse Trendelenburg

Bring beauty and character into your room with exceptional bed functions, made in Australia.

AMSTERDAM

Headboard and footboard























AFTER DARK FEATURES



Flash light on the back of the remote



Under bed LED lighting



Two personal programmable memory position settings



Zero gravity preset functions



Patient transfer pre-set function



Backlit remote



- Backlit remote activated by a touch of a button.
- Remote preprogrammed functions:
 - M1 user personal settings.
 - M2 user personal settings.
 - ZG Zero Gravity preset functionoverride for personal settings.
 - PT Patient Transfer preset functionoverride for personal settings.
- · LED under bed lighting.
- Flashlight on/off button.
- Indicator light: red- battery low, greencharged

BED CODES & SIZES

Sizes for every room size and lifestyle.

BED SIZE	CODE	MATTRESS SURFACE SIZE	BASE SIZE
Long Single	IC777LS	2030 x 900mm	2150 x 960mm
King Single	IC777KS	2030 x 1070mm	2150 x 1120mm
Long Double One hand control	IC777D	2030 x 1350mm	2150 x 1400mm
One Piece Queen	IC777Q	2030 x 1500mm	2150 x 1580mm
Split Queen	IC777SQ	2030 x 750mm (per side)	2150 x 1580mm
Split King	IC777K	2030 x 900mm (per side)	2150 x 940mm

NOTES:

SPLIT QUEEN

- · Separate head and knee functions
- · One piece hi-lo and trendelenberg
- · Individual hand control per side
- · Fits a one piece queen headboard

SPLIT KING

Always remember that TWO Long Single beds together make a King size bed.

TECHNICAL SPECS

Key details and information.

- 200kg safe working load on long single and king single sizes and 300kg for double and queen sizes. The safe working load includes any model of Icare mattresses.
- Long-life rechargeable battery back up, lasting 2-5 days of usual use.
- Easy press button remote, backlit when activated for use.
 Wipeable.
- Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Strong metal mesh mattress platform allowing increased airflow
- Actuators 1x 6000 newton 24 watt head adjust actuator. 1x 6000 newton - 24 watt knee adjust actuator. 2x 6000 newton - 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- No pull motors for added safety.
- Wheels 8x heavy duty swivel castors with one lockable castor on each corner.



HEADBOARDS & FOOTBOARDS

Codes and measurements.





AMSTERDAM HEADBOARD

LONG SINGLE Code: HB7AS Width: 960mm

KING SINGLE Code: HB7AKS Width: 1120mm

DOUBLE Code: HB7AD Width: 1400mm

QUEEN Code: HB7AQ Width: 1580mm

COLOGNE HEADBOARD

LONG SINGLE Code: HB7CS Width: 960mm

KING SINGLE Code: HB7CKS Width: 1120mm

DOUBLE Code: HB7CD Width: 1400mm

QUEEN Code: HB7CQ Width: 1580mm

KARLSTEAD HEADBOARD

LONG SINGLE Code: HB7KS Width: 960mm

KING SINGLE Code: HB7KKS Width: 1120mm

DOUBLE Code: HB7KD Width: 1400mm

QUEEN Code: HB7KQ Width: 1580mm

AMSTERDAM FOOTBOARD

LONG SINGLE Code: FB7AS Width: 960mm

KING SINGLE Code: FB7AKS Width: 1120mm

DOUBLE Code: FB7AD Width: 1400mm

QUEEN Code: FB7AQ Width: 1580mm

COLOGNE FOOTBOARD

LONG SINGLE Code: FB7CS Width: 960mm

KING SINGLE Code: FB7CKS Width: 1120mm

DOUBLE Code: FB7CD Width: 1400mm

QUEEN Code: FB7CQ Width: 1580mm

KARLSTEAD FOOTBOARD

LONG SINGLE Code: FB7KS Width: 960mm

KING SINGLE Code: FB7KKS Width: 1120mm

DOUBLE Code: FB7KD Width: 1400mm

QUEEN Code: FB7KQ Width: 1580mm

ICARE IC777

A NEW WAY TO CARE





ICARE IC777 PREMIUM HOMECARE BED





ICARE MEDICAL GROUP

icaremedicalgroup.com.au info@sleepsystems.com.au

1800 105 853 | 03 9070 2220