



NEW PRODUCT

HOIST MAT

The Hoist mat surface is ultra clear and made from twice-extruded heavy duty PVC, designed for mechanical hoists, walkers and heavy load bearing. The Hoist mat is secured to the floor with rounded grippers or an adhesive backing.

The Hoist mat creates a strong flat stable surface on top of a carpeted floor that enables you to place a mechanical hoist on top and to safely lift patients in and out of hospital beds.



IDEAL FOR HOSPITALS, AGED CARE AND PATIENTS PRIVATE ROOMS

Product Construction	Twice-extruded heavy duty PVC
Surface	Non-slip orange peel texture
Standard Edging	Tapered anti trip edges – all sides
Cleaning	Mop with water, soap, mild detergents or vacuum
Usage Areas	On carpet - patients private rooms / hospital suites / aged care
Mat Size	1200mm x 2400mm - Rectangular
Mat Profile	3.5mm
Colours	Ultra Clear

rehab**hire**&sales

310 Lorimer Street Port Melbourne Victoria 3207
t: 03 9646 7200 | e: contact@rehabhire.com.au | www.rehabhire.com.au