



SMOOV

by Alber GmbH

MADE IN GERMANY

max speed	10 km/h
user capacity	140kg
range	up to 20km
battery	6.2Ah, 36v, 224Wh
motor	250W / 450W peak 12Nm / 25Nm peak
unit weight	7.2kg

- swiveling drive wheel enables all terrain use
- rigid and folding frame compatibility
- stop with mounted button
- stop with pushrims at low indoor settings
- adjust speed without disengaging
- easy attachment and removal
- USB C output to charge devices from smooov
- connect to smart phone with free app
- safe for air travel according to IATA

www.smooovbyalber.com.au

