

NEW


IC30 ACTIVEX™

icare
Medical Group

The IC30 Gel-Infused ActiveX™ Mattress with Charcoal and Firm Edge Support is designed to provide comfort and support for a good night's sleep. The gel infusion in the mattress helps to regulate temperature, ensuring a cool and comfortable sleep surface. The addition of charcoal helps to absorb odors and moisture, promoting a clean and hygienic sleeping environment.

ActiveX™ is temperature and pressure sensitive which provides an optimum balance of immersion and envelopment creating a sleep surface with excellent support, comfort and very high level of pressure care.

The firm edge support helps to maintain the shape and structure of the mattress, reducing the risk of roll-out and assisting with transfer access. This also makes it ideal for use on profiling beds, which are beds that can be adjusted to various positions to provide comfort and support for those with limited mobility.



SPECIFICATIONS

Feel Rating	Medium
Mattress Height	200mm
Risk Category	High Risk
User Weight	20kg - 228kg
Warranty	10 Year Warranty

BED CODES & SIZES

LONG SINGLE	IC30LS	2030 X 900mm
KING SINGLE	IC30KS	2030 x 1070mm
LONG DOUBLE	IC30LD	2030 x 1350mm
QUEEN	IC30Q	2030 x 1520mm
SPLIT QUEEN	IC30SQ	2030 x 760mm (per side)

PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <16.5mmHg.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For more information refer to the Icare pressure care manual.