

# Deluxe-Flex

Gel-Infused Memory Foam Mattress

**Premium*Lift***




# Enjoy a healthier, more comfortable sleep environment

## Deluxe-Flex™

Designed to improve your sleep...

Sleep can change as we age, it is important to emphasize that there are many benefits of getting a good night's sleep.

Older adults who regularly get a good night's sleep report improved cognitive function, improved memory, and better mental health. This is important for preventing age-related cognitive decline and dementia, as well as play a role in preventing depression, which is a common mental health disorder among older individuals.



**TGA** Approved  
Medical Device

and keep you healthy

Using an innovative combination of Ice fiber breathable cover and gel-infused memory foam, Deluxe-Flex™ will regulate your body temperature and balance your body weight across the mattress uniformly.

Chronic poor sleep can significantly affect our bodies and minds. Studies have shown that poor sleep habits can increase the risk of dementia by 33%. It can also increase risks of heart disease, type 2 diabetes, high blood pressure, obesity, and age our brain 3-5 years. Most importantly, it can decrease our immunity which makes it hard for our body to fight off infections. We also are three times more likely to catch a cold. So, let's fix our sleep and keep our body healthy.

**PremiumLift**

Advanced quality you can trust

**PremiumLift**



## Deluxe-Flex™

### Technical Information



#### Clinical Features

- ✓ Gel infused top layer
- ✓ High density foam base
- ✓ Firm walls
- ✓ Ice fiber breathable & washable cover
- ✓ Fire retardant foam
- ✓ 360° / two-sided integrated zipper
- ✓ Inner zipped washable cover will provide extra hygiene
- ✓ Superior durability
- ✓ Excellent performance

#### Compliance

ARTG 398711

ISO13485

Australian Standards

BS7177, BS5852



Firmer sides improve comfort and support while entering and exiting bed



Gel Memory Foam regulate your body temperature and improve your sleep

#### Product Code

PLDFKS (King Single)

#### Risk Factor

Low Risk

#### Maximum Patient Weight

220 kg

#### Mattress Dimensions

King Single: 203cm x 106cm x 20cm

#### Mattress Weight

15.50 kgs

#### Colour

White top cover

Grey bottom cover

#### Applications

Homecare mattress

#### Warranty:

2-year cover

5-year foam