


Multi-Flex

Advanced Pressure Care Mattress

Premium*Lift*



Multi-Flex is suitable for every care environment



Multi-Flex Superior Pressure Care Reduction

The clinically proven Multi-Flex™ pressure reducing mattress has been developed using the very latest manufacturing technology. This high quality product features a Gel-infused castellated foam surface, ensuring that total comfort is achieved whilst delivering excellent pressure re-distribution.

The pressure reducing, gel memory foam castellated zones are designed to maximise pressure re-distribution in the most vulnerable areas, regulate body temperature and improve sleep quality.

The castellations not only provide superior pressure relief and postural support but the channels cuts ensure that the mattress provides a constant airflow around the patient. This airflow assists in controlling temperature and therefore endeavours to improve comfort and the patients ability to sleep.

The high-density foam base and walls offers total stability during side transfer procedure, to alternative seating systems or that of assisting stand transfer.

TGA Approved
Medical Device

PremiumLift


High Quality, Superior Durability, Medical Grade, Pressure Care, Pain Relief

Multi-Flex™


Features & Benefits



4-Way stretch, vapour permeable cover protected by a cotton inner lining



Gel-infused castellated foam surface, ensuring that total comfort is achieved whilst delivering excellent pressure re-distribution



Clinically proven pressure reducing mattress developed using the very latest manufacturing technology



High density foam base and firm hinged walls offers total stability during side transfer



Waterfall fluid penetration flap and hardwearing base cover protects foam core against the bed platform



Gel memory foam castellated top layer regulates body temperature and improve sleep quality

PremiumLift

Advanced quality you can trust

PremiumLift



Multi-Flex™

Technical Information



Clinical Features

- ✓ Castellated gel infused top layer
- ✓ High density foam base
- ✓ Hinged firm walls
- ✓ 4-Way stretch antimicrobial cover
- ✓ Waterproof and resistant to patient body fluids
- ✓ Resistance to cleaning agents
- ✓ Wipe clean
- ✓ Stain resistant
- ✓ Heavy duty PVC base
- ✓ Fire retardant foam
- ✓ 360° / two-sided integrated zipper
- ✓ High frequency sealed & welded seams
- ✓ Waterfall fluid penetration flap
- ✓ Superior durability
- ✓ Excellent performance

Compliance

ARTG 398711

ISO13485

Australian Standards

BS7177, BS5852

Warranty:

2-year cover

5-year foam

Product Code

PLMFS (Single)

PLMFKS (King Single)

Risk Factor

Medium / High Risk

Maximum Patient Weight

220 kg

Mattress Dimensions

Single: 198cm x 88cm x 15cm

King Single: 198cm x 106cm x 15cm

Mattress Weight

13.98 kgs

Colour

Purple top cover

Grey bottom cover

Applications

General medical mattress

Aged care mattress

General hospital ward mattress

Homecare mattress

premiumlift.com.au