



VIDA WHEELCHAIR

The Aspire VIDA is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends. Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support. The 21 stone safe working load caters to a broad range of users.



UPHOLSTERY
Orange/Black
FRAME
Silver



UPHOLSTERY
Green/Black
FRAME
Silver



5 YEARS WARRANTY
Frame

1 YEAR WARRANTY
Components

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk

FEATURES & BENEFITS

- Available as attendant propelled or self-propelled
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Rear storage pocket for essentials
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)

Self - Propelled



Attendant Propelled



SWL
135kg
21st

	OVERALL LENGTH	OVERALL WIDTH	OVERALL HEIGHT	SEAT DEPTH	SEAT WIDTH	ANTERIOR SEAT HEIGHT	POSTERIOR SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	REAR WHEEL	MIN. CHAIR WEIGHT	TOTAL CHAIR WEIGHT
	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	in	kg/lbs	kg/lbs
Attendant Propelled	995 / 39	655 / 26	930 / 36.5	420 / 16.5	450 / 18	495 / 19.5	445 / 17.5	430 / 17	200 / 8	13	8.8 / 19.4	9.6 / 21.2
Self Propelled	1030 / 40.5	685 / 27	930 / 36.5	420 / 16.5	450 / 18	495 / 19.5	445 / 17.5	430 / 17	200 / 8	22	7.5 / 16.6	11.7 / 25.8

ASFY99,7 141123