



Barwon

Hudson

Order Form: Barwon & Hudson Lift Chair

Prepared by

Client Name

Order Date

Reference No.

Circle required dimensions and configurations

CONFIGURATION

Mechanism	Single Motor (Standard)	Dual Motor	180kg Dual Motor RC	Electric Recline RC	Manual Recline CP	Rocker Manual Recline CP
Wheels	50mm Castors (Standard)	75 Mobile Castors (+25mm Seat height)			Glides	Remove Castors on Site
Memory Foam	Seat	Legrest		Backrest		Armrest
Back Cushion Options	Gluteal Recess	Concave Back		Blow-fill Lumbar		Comfort Pack
ROHO Cushion Insert	None (Standard)	18x16 Low Profile	18x16 High Profile	16x16 Low Profile	16x16 High Profile	Include Infill Plug
Controller Option	Right (Standard)	Left		Double Pocket		Joystick
Back Angle	Standard		Vertical			Over-Vertical
Other Options	Waterfall Arms		Scotch Guard			Tablet Arm Left/Right
	Flush Front		Leg rest Side Walls			Controller Loop Left/Right
	Full chair cover		Velcro armrest/headrest covers		Loose armrest/headrest covers	
	Removable Headrest A B C D E (please circle one)			24V Full Battery Backup		

DIMENSIONS

	Seat Depth (mm)	Seat Height (mm)	Back Height (mm)	Footrest Length (mm)	Scapular Height (mm)	Seat Width (mm)	Arm Height (mm)	Seat Foam (mm)
Decrease Size	-100 CP -150 CP -25 -50 -75	-30 RC -25	-120 CP -150 CP -180 CP -30 CP -60 CP -90 CP	-50 -25		-150 -50 -100	-60 -30	
Mini Standard Size	455	456	636	N/A	429	540	199	N/A
Increase Size	+25 +50 +75	+20	+30 +60 +90 +120 +150 DM	+25 +50		+50	+30 +60	+20
Decrease Size	-100 -150 CP -25 -50 -75	-30 RC -25	-120 CP -150 CP -180 CP -30 -60 -90 CP	-50 -25		-150 -50 -100	-60 -30	-20
A Standard Size	490	476	647	N/A	440	540	180	N/A
Increase Size	+25 +50 +75	+20	+30 +60 +90 +120 DM +150 DM	+25 +50		+50	+30 +60	+20
Decrease Size	-100 -150 CP -25 -50 -75	-30 RC -25	-120 CP -150 CP -180 CP -30 -60 -90 CP	-50 -25		-150 -50 -100	-60 -30	-20
B Standard Size	520	496	688	N/A	481	540	190	N/A
Increase Size	+25 +50	+25	+30 +60 DM +90 DM +120 DM +150 DM	+25 +50		+50	+30 +60	+20
Decrease Size	-100 -150 -25 -50 -75	-30 RC -25	-120 CP -150 CP -180 CP -30 -60 -90	-50 -25		-150 -50 -100	-60 -30	-20
C Standard Size	545	546	718	N/A	511	540	190	N/A
Increase Size	+25	+25 +50	+30 DM +60 DM +90 DM +120 DM	+25 +50		+50	+30 +60	+20

NOTE: Reduced seat foam limits Safe working load to 90kg and reduces seat foam warranty to 2 years.

* DM = Dual Motor Only * RC = Remove Castors -30mm Seat Height * CP = Not compatible with Comfort Pack

COLOUR AND FINISH

Fabric	Inner					
	Outer					
Show wood		Walnut	Mahogany	Baltic	Teak	Remove Show Wood

ADDITIONAL NOTES

