



aspire

LIFESTYLE COMMUNITY BEDS

Adjustable and stylish Homecare
beds with a modern look suited to the
Homecare Environment.



www.aspirecare.com.au
1300 133 120



LIFESTYLE COMMUNITY BED



The Aspire Lifestyle Community Bed range has been designed to seamlessly blend into contemporary Homecare environments.

- Available in Single & King Single sizes
- Tool free folding and unfolding feature
- Electronics remain in place during assembly/disassembly, folding and transportation
- Easy handling and movement on castors with transport brackets included
- Built to withstand the rigors of regular transportation & assembly

HIGHLY ADJUSTABLE


Four electric profiling functions to suit all comfort and positioning needs.

BED EXTENSION




150mm

BED HEIGHT RANGE




245 - 645mm

BACKREST INCLINE/DECLINE




70°

KNEE BEND INCLINE/DECLINE



30°

TRENDELENBURG TILT



10.5°



Pictured with optional bed end shrouds and full length side rails

STYLISH TIMBER END SHROUDS

Sophisticated timber panelling to suit the homecare environment.

Optional Accessory

Single [BEA029850](#)
King Single [BEA029855](#)

FOLDING SIDE RAILS

3/4 length fold down side rails ensure user safety and fold conveniently out of the way as desired.

Optional Accessory

[BEA016586](#)

DURABLE LAMINATE SIDE RAILS

Sleek timber panelled side rails with multiple positioning heights can be lowered out of the way as desired.

Optional Accessory

Side Rails Full Length [BEA029863](#)
Side Rails Extension Kit [BEA029864](#)

LEADING QUALITY FRAME & ELECTRONICS

Strong durable steel frame with Linak electronics and handset allows for market leading warranty provision.

FOUR TOTAL-LOCK CASTORS

Safely locate the bed and prevent castor roll and directional rotation.

TRANSPORTABLE DESIGN

Tool free with no requirement to remove electronics during assembly and set-up, ensuring quick and easy readiness upon delivery or for transportation and storage.

Bed comes pre-assembled in the upright transport configuration shown, including transport brackets. Simply wheel the bed from the carton.

BEBO46030

Single

BEBO46040

King Single



Safe Working Load

OTHER OPTIONAL ACCESSORIES



Self Help Pole
[BEA016596](#)



IV Pole
[BEA016592](#)

MAXI LIFESTYLE COMMUNITY BED




The Aspire Maxi Lifestyle Community Bed is a stylish, convenient and economical solution for users weighing up to 350kg.

- Unique power assisted electric folding and opening feature
- Safe moving and handling for carers, staff and distributors
- Highly maneuverable
- Small storage footprint
- Easy transportation
- Holistic solution for bariatric care

HIGHLY ADJUSTABLE


Four electric profiling functions to suit all comfort and positioning needs.

BED HEIGHT RANGE



330 - 730mm

BACKREST INCLINE/DECLINE



70°

LEG RAISE INCLINE/DECLINE



30°

TRENDELENBURG TILT



10°



STYLISH TIMBER BED ENDS


Offer a simple solution for cleaning and infection control.



REMOVABLE FOLDING SIDE RAILS

Durable design to handle the rigors of equipment rental pools.

Optional Accessory
BEA016586



SUPER-KING SIZE MATTRESS PLATFORM

Reduces falls risk and moving and handling complexity.

Compatible with the Aspire BARI 10 alternating turning mattress offering a high care rental package at a market leading price.



AUTO-REGRESSION IN THE BACKREST

Reduces skin shear forces and removes a significant cause of patient slippage. This reduces the need for care staff to reposition oversized patients.




LEADING QUALITY FRAME & ELECTRONICS

Linak electronics offer market leading reliability and performance.



RAPID ELECTRIC FOLDING

Single button press electrical folding mechanism requires no lifting and improves workplace health & safety outcomes for delivery staff and maintenance teams.



LARGE 120MM CASTORS

Improve ease of movement during transportation and storage.

OTHER OPTIONAL ACCESSORIES



Maxi Self Help Pole
BEA016598



IV Pole
BEA016592

BEB046020

MAXI - Super King



Safe Working Load

SPECIFICATIONS



Size	Single	King Single	MAXI - Super King
Code	BEB046030	BEB046040	BEB046020
Height Range	245 - 645 mm		330 - 730 mm
Overall Length	2220 - 2370 mm		2290 mm
Overall Width	995 mm	1155 mm	1360 mm
Recommended Mattress Size (L x W x H)	2000 x 900 x 150 mm	2030 x 1050 x 150 mm	2060 x 1150 x 200 mm
Folded Dimensions (L x W x H)	405 x 995 x 1165 mm	405 x 1155 x 1165 mm	415 x 1360 x 1270 mm
Backrest Incline	70°		73°
Knee Bend	30°		30°
Maximum Patient Weight	180 kg		350 kg
S.W.L	220 kg		385 kg
Product Weight	75.5 kg	83 kg	155 kg
Compliances	ARTG AUST R ARTG 239178, IEC60601-2.52:2010		ARTG AUST R ARTG 239178, IEC60601-2.52:2010


1 YEAR WARRANTY
Handset

2 YEAR WARRANTY
Electronics

3 YEAR WARRANTY
Frame

FOR HOMECARE OR RENTAL USE

Aspire Lifestyle Community Beds can be dressed up with bed end shrouds and side rails for the homecare environment or be provided as standard for the rental market.



SLEEP SURFACES

OUR RANGES


LIFECOMFORT® Foam Mattresses

All **LIFECOMFORT®** Mattresses are fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments. Multiple foam layering combinations provide maximum comfort, support and pressure redistribution across a wide range of user shapes and sizes.

HYBRID® & ACTIVE AIR® Dynamic Alternating Air Mattresses

Aspire alternating air mattresses are designed to offer comfort and support for patients at high risk of pressure sores.

FOR THE LIFESTYLE COMMUNITY SINGLE AND KING SINGLE BEDS




**Lifecomfort
Standard**

BEM046410 - Single
BEM046460 - King Single




**Lifecomfort
Deluxe**

BEM046510 - Single
BEM046455 - King Single




**Lifecomfort
Pressure Reducing**

BEM046610 - Single
BEM046660 - King Single



**Lifecomfort
Premier Care**

BEM046725 - Single
BEM046730 - King Single



Hybrid Comfort

PTM638005 - Single
PTM638006 - King Single



Active Air 8

PTM639300 - Single
PTM639310 - King Single

FOR THE MAXI LIFESTYLE COMMUNITY BED



**Lifecomfort
Maxi Lifestyle**

BEM046710 - Super King



Active Air BARI 10

PTM639520 - Super King

ASSOCIATED PRODUCTS



MOBILITY SCOOTERS

Wide range of Aspire Scooters for all user needs.

Shown Aspire HS268 Manual Folding Mobility Scooter SCT671135



LIFT CHAIR RECLINERS

Luxury Aspire recliner range has four, dual and single motor functions.

Shown Aspire Raphael Quattro Powerlift Recline Chair Sizes S-XL Available



SEAT WALKERS

Our range includes various options to suit different conditions, abilities and personal requirements.

Shown Aspire Vogue Carbon Fibre Seat Walker WAF705350



OVERBED TABLES

Lightweight mobile table positions easily over the lifestyle community beds.

Shown Over Bed Table with Gunmetal Frame & Sublime Teak Laminate Top BEA006705