



via GO[®]

The compact push aid for wheelchairs

NEW



Pushing
the easy way.

via GO – push smart, not hard.



The powerful, easy-to-use push aid for wheelchairs.

There are times when even a poor pavement surface or a slight incline can become an insurmountable obstacle for everyday tasks. via GO – the powerful, easy-to-use push aid from Alber – always provides the right amount of assistance for the accompanying person, helping to effortlessly master such situations with confidence. Rediscover independent mobility and enjoy any excursion in a relaxed manner.

The high-quality and particularly compact push aid fits almost all common wheelchair models. Your trusted medical supply dealer will check whether your wheelchair is compatible and provide expert support for fitting.

via GO – pushing the limits.

COMPACT AND CONVENIENT



Experience more:
simply practical.

With its impressive versatility, via GO will become your everyday companion – for day-to-day tasks and excursions alike. Thanks to its compact dimensions and lightweight design, the push aid is very easy to handle and can be stowed away quickly to save space.

STRONG AND ENDURING



No effort:
simply move on.

Every town has inclines and ramps. The via GO gives you the right amount of assistance just when you need it, helping you to get your loved ones on the move with no problem at all – for up to 10 miles.

INTUITIVE AND FLEXIBLE



Clever thinking:
simply done.

Even without any previous knowledge, the via GO is so straightforward you can use it immediately. Large buttons with a clear design make control a breeze. What's more, you don't have to remove the via GO to fold up your wheelchair.

SAFE AND RELIABLE



Alber quality:
simply high class.

The via GO push aid incorporates the depth and breadth of experience and expertise of Alber, the specialist and market leader in add-on drives for wheelchairs. Alber products regularly set new benchmarks and stand for enduring reliability and TÜV-tested safety.



MADE IN GERMANY 

Technical data for via GO

Dimensions

Drive unit: 6.7 x 9.4 x 10.6 in (W x H x L)

Control unit: 1.4 x 3.7 x 1.1 in (W x H x D)

Speed

0.9–3.4 mph

Range

up to 10 miles

integrated rechargeable

lithium-ion battery*

Weight capacity

up to 330 lbs (user weight)**

Gradient

up to 16%**

Weight

Charger: 1.19 lbs

Drive unit: 12.87 lbs

Control unit: 0.31 lbs

* Depending on user weight, topography and ambient temperature in accordance with ISO 7176-4

** May be lower depending on the wheelchair used and operator ability




Daily life: simply easier.

Both the control and the installation of the components are particularly uncomplicated and user-friendly with via GO.

The small, ergonomic **control unit** attaches simply to the push handles of your wheelchair.

The clever **mounting device** allows your wheelchair to be folded up as normal, even with the push aid still attached.

As a result of its lightweight design and large handle, the **drive unit** is extremely easy to attach to and remove from the wheelchair.



www.alber.de

via GO[®]

The compact push aid for wheelchairs



28.0004.6.02.01 Version 03/2023 © Alber GmbH, Albstadt.
Subject to technical changes.

At alber.de you can find even more information about via GO or about reimbursement by your health insurance company. The customer service team at Alber GmbH looks forward to meeting you and will be happy to provide you with advice and assistance.

The via GO is CE compliant according to the European Medical Device Regulation (MDR) 2017/745 and has been tested on a voluntary basis by TÜV Süd according to the latest edition of the product standard EN 12184.

Our quality commitment: your safety! Alber is certified under the current quality management standard ISO 13485 for medical devices. Our environmental management system, which is certified under ISO 14001, ensures the environmentally responsible manufacture of our products.

Alber GmbH

Vor dem Weißen Stein 14
72461 Albstadt-Tailfingen, Germany
info@alber.de

Phone: +49 (0)7432 2006-0

www.alber.de



LET'S via GO!

rehabhire & sales

320 Lorimer Street Port Melbourne Victoria 3207
t: 03 9646 7200 | e: contact@rehabhire.com.au | rehabhire.com.au

MADE IN GERMANY 