

Heavy Duty Rollator

Features:

- ✓ Lightweight aluminium frame with soft touch brakes
- ✓ Durable 8" castors
- ✓ Easy to fold with removable backrest for storage and transportation
- ✓ Padded seat with storage bag under seat

**Tech Specs**

Seat width	44 cm
Seat depth	39 cm
Seat height	54-56 cm
Overall width	67 cm
Overall depth	81 cm
Handle height	84-94 cm
Weight capacity (static)	180 kg
Total weight	9.5 kg