



FLEX SEAT WALKER



6" Wheels

8" Wheels



Champagne



Black



Lilac



Blue

- Super lightweight making it easy to manoeuvre and easy to transport
- Available in 6" and 8" wheel variations
- Robust aluminium frame for improved durability and stability
- Height adjustable seat height to suit most users
- Soft padded seat for improved seating comfort
- Adjustable and lockable brake handles



FLEX	Overall Length mm/ in	Overall Width mm/ in	Overall Height mm/ in	Seat Depth mm/ in	Seat Width mm/ in	Seat Height mm/ in	Backrest Height mm/ in	Folded Width mm/ in	Product Weight kg/ lbs	Castor Size in
6" Wheels WAF709660CH Champagne WAF709660PU Lilac	730 / 28.7	570 / 22.4	730 - 830 / 28.7 - 32.6	300 / 11.8	360 / 14.1	430 - 520 / 16.9 - 20.4	710 / 27.9	290 / 11.4	6.3 / 13.8	6
8" Wheels WAF709670BK Black WAF709670BL Blue	820 / 32.2	585 / 23	780 - 970 / 30.7 - 38.1	300 / 11.8	360 / 14.1	475 - 570 / 18.7 - 22.4	760 / 29.2	390 / 15.3	7.3 / 16	8

ASFY9715 031224